

Dialysis Tests

DiaQuant[®]
Test strips for dialysis



Decades of experience in test strips

✓ **Reliable and accurate**

Designed for dialysis centres

✓ **Convenient and effective**

Dip & read procedure

✓ **Simple and fast**

MACHEREY-NAGEL
www.mn-net.com



DiaQuant® Test strips for dialysis

■ DiaQuant® Peracetic Acid

Test for residual peracetic acid in rinse water



This test is specific to peracetic acid. The easy Dip & Read procedure provides reliable results in 10 seconds.

Compared to peroxides, peracetic acid is the more powerful disinfectant. It is often used in the disinfection of dialysate lines and dialysers. Careful rinsing of the system is mandatory, before the next patient is connected. Both, DiaQuant® Peracetic acid and DiaQuant® Peroxide may be used to ensure, that the system is free of disinfectant residues. DiaQuant® Peracetic acid, however, is specific to peracetic acid

REF: 932 001
Type: test strip
Range: 0 – 5 – 10 – 20 – 30 – 50 ppm (mg/l) peracetic acid (PAA)
Sufficient for: 100 tests
Shelf life: at least 2.5 years after production
Storage: 4 – 30 °C (39 – 86 °F)
Accessories: no accessories are necessary to perform this test



■ DiaQuant® Peroxide

Test for residual peroxide in rinse water



This test is for the quick and easy detection of Peroxide in rinse water. The easy Dip & Read procedure provides reliable results in 10 seconds.

Health care professionals in dialysis centers regularly clean dialysate lines or dialyzers with sanitizers based on peroxides or peracetic acid. They need to ensure, that the disinfectant is completely rinsed from the system prior to the next use.

DiaQuant® Peroxide provides a quick and easy means to detect peroxides in rinse water down to 1 ppm (mg/l) in order to verify, that no disinfectant residues remain in the system.

REF: 932 002
Type: test strip
Range: 0 – 1 – 3 – 10 – 30 – 100 ppm (mg/l) H₂O₂
Sufficient for: 100 tests
Shelf life: at least 2.5 years after production
Storage: 4 – 30 °C (39 – 86 °F)
Accessories: no accessories are necessary to perform this test



DiaQuant® Test strips for dialysis

■ DiaQuant® pH 3.6 – 6.1

Test for rinse water after hot disinfection with hydroxy acetic acid containing disinfectants




Depending on the disinfection demands, health care professionals use sanitizers based on hydroxy acetic acid. Before the next patient is admitted, the system has to be checked for sanitizer residues. For hydroxy acetic acid based disinfection, DiaQuant® pH is the ideal control. The easy Dip & Read procedure provides dependable results in 5 seconds. When the pH-strip does not detect any acid, the system is completely rinsed.

REF: 932 003
Type: test strip
Range: pH 3.6 – 4.1 – 4.4 – 4.7 – 5.0 – 5.3 – 5.6 – 6.1
Sufficient for: 100 tests
Shelf life: at least 2.5 years after production
Storage: 4 – 30 °C (39 – 86 °F)
Accessories: no accessories are necessary to perform this test



■ DiaQuant® Chlorine

Sensitive test for total chlorine including chloramines



DiaQuant® Chlorine allows the quick and easy determination of chlorine in water. The easy Dip & Read procedure provides reliable results within 10 seconds.

When the status of dialysis machines or dialysate lines indicates the use of powerful disinfecting agents, health care professional often use chlorine-based disinfectants. Before the next use, they need to check the system carefully for residual chlorine. DiaQuant® Chlorine provides safe and easy means to detect chlorine residues.

Low chlorine levels are also important for feed water, which is purified using reverse osmosis (RO) systems. Chlorine affects the RO unit and may lead to early breakdown of the membrane. DiaQuant® Chlorine allows the safe and easy detection of chlorine including chloramines in water down to 0.1 ppm (mg/l). It helps to ensure the longevity of expensive RO units.

REF: 932 004
Type: test strip
Range: 0 – 0.1 – 0.5 – 1 – 3 – 10 ppm (mg/l) Cl₂
Sufficient for: 100 tests
Shelf life: at least 2.5 years after production
Storage: 4 – 30 °C (39 – 86 °F)
Accessories: no accessories are necessary to perform this test



DiaQuant® Test strips for dialysis

DiaQuant® Water Hardness

Test for low level (5 ppm) water hardness in feed water



This test allows the convenient and reliable determination of low-range water hardness in feed water. The easy Dip & Read procedure provides reliable results in 10 seconds.

DiaQuant® Water Hardness is used after the water-softener stage of the clean-up process. Softening is a pretreatment for water that usually undergoes further purification processes like reverse osmosis. Water hardness is the total concentration of calcium and magnesium expressed in degree of water hardness or in mg/l (ppm). With DiaQuant® Water Hardness levels down to 0.3°d or 5 ppm (mg/l) can safely be detected ensuring proper feed water for the reverse osmosis system.

- REF: 932 005
Type: test strip
Range: 0 – 5 – 10 – 20 ppm (mg/l) CaCO₃
Sufficient for: 100 tests
Shelf life: at least 2.5 years after production
Storage: 4 – 30 °C (39 – 86 °F)
Accessories: no accessories are necessary to perform this test



DiaQuant® Peracetic Acid P

Test for potency levels of peracetic acid in dialyzer reprocessing.



This test is specific to peracetic acid and it will not detect peroxides. The easy Dip & Read procedure provides reliable results in 10 seconds.

The test is a convenient and accurate way to check the concentration of peracetic acid in dialyzer reprocessing. It detects potency levels up to 2000 ppm (mg/l) without the need for dilution, which makes the testing safe and easy. DiaQuant® Peracetic acid is specific for peracetic acid and does not react with peroxides.

- REF: 932 006
Type: test strip
Range: 0 – 500 – 1000 – 1500 – 2000 ppm (mg/l) peracetic acid (PAA)
Sufficient for: 100 tests
Shelf life: at least 2.5 years after production
Storage: 4 – 30 °C (39 – 86 °F)
Accessories: no accessories are necessary to perform this test



DiaQuant® Test strips for dialysis

DiaQuant® test strips

Dip & Read tests for dialysis centres







In dialysis, reliable on-site chemical testing is a critical factor for patient safety. Accurate and quick testing is also fundamental to run a dialysis center effectively and successfully. The new DiaQuant® line from MACHEREY-NAGEL is designed for the dialysis environment. Our proven technology makes testing reliable and accurate, thus providing a safe environment for your patients.


Fast results
Health care professionals use DiaQuant® tests for a variety of different substances. The product group offers a complete range of testing equipment for disinfection and water supply. All tests follow an easy dip&read procedure and provide reliable results in 5 – 10 seconds.

Easy to use
All DiaQuant® tests are ready-to-use kits. They are pre-calibrated and do not need further equipment or reagents. As “lab in the pocket” they are easy to use for any health care professional.

Precise readings
The colour charts are adjusted and checked using standard solutions. Users can be sure to receive accurate readings whenever they test.



Test	Range	Rinse Water	Feed Water	Potency Testing	REF
 DiaQuant® Peracetic Acid	0 – 5 – 10 – 20 – 30 – 50 ppm (mg/l) PAA	✓			932001
 DiaQuant® Peroxide	0 – 1 – 3 – 10 – 30 – 100 ppm (mg/l) H ₂ O ₂	✓			932002
 DiaQuant® pH 3.6 – 6.1	pH 3.6 – 4.1 – 4.4 – 4.7 – 5.0 – 5.3 – 5.6 – 6.1	✓			932003
 DiaQuant® Chlorine	0 – 0.1 – 0.5 – 1 – 3 – 10 ppm (mg/l) Cl ₂ (total)	✓	✓		932004
 DiaQuant® Water Hardness	0 – 5 – 10 – 20 ppm (mg/l) CaCO ₃		✓		932005
 DiaQuant® Peracetic Acid P	0 – 500 – 1000 – 1500 – 2000 ppm (mg/l) PAA			✓	932006

 All DiaQuant® Test strips are CE-marked according to the medical device directive 93/42/EWG

DiaQuant® Test strips for dialysis

■ MACHEREY-NAGEL - Your competent partner in analytical chemistry

Tradition and Modernity – over 90 years experience

- Worldwide operating German company (founded 1911)
- Subsidiaries in France, Switzerland and USA
- Distributors in over 150 countries
- Longtime tradition in filter papers
- Development, production and sales of special products for water, environmental and food analysis, for biotechnology, chemical and pharmaceutical industry and medical diagnostics

Quality and Diversity – 5 product ranges with over 25 000 products „Made in Germany“

Filtration	Filter Papers Filter Membranes Extraction Thimbles
Rapid Tests	Test Papers and Test Strips Urine Test Strips
Water Analysis	Colorimetric and Titrimetric Test Kits Photometric Water Analysis
Chromatography	HPLC GC TLC SPE
Bioanalysis	Kits for Purification of Nucleic Acids and Proteins Transfer Membranes

Service

- Technical support by our Customer Service Center
- Product-specific quality certificates and Material Safety Data Sheets
- Customised special products
- Easy ordering and short delivery times
- Large number of catalogues, flyer and product information – of course free of charge

Contact

Technical Support and Customer Service Center for
Water Analysis, Filtration, Rapid Tests, Medi-Test:

+49 (0) 24 21 969 161

+49 (0) 24 21 969 174

Product Management Filtration and Rapid Tests:

+49 (0) 24 21 969 166

Please visit our website:

www.mn-net.com

Your local distributor

www.mn-net.com

KATEN100103/FlyerDiaQuant en1/5/0/11.2009 PD
Printed in Germany

MACHEREY-NAGEL

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6 – 8 · 52355 Düren · Germany

Germany
and international

Tel.: +49 (0) 24 21 96 90

Fax: +49 (0) 24 21 96 91 19

e-mail: sales-de@mn-net.com

Switzerland
MACHEREY-NAGEL AG

Tel.: +41 (0) 62 388 55 00

Fax: +41 (0) 62 388 55 05

e-mail: sales-ch@mn-net.com

France
MACHEREY-NAGEL EURL

Tel.: +33 (0) 3 88 68 22 68

Fax: +33 (0) 3 88 51 76 88

e-mail: sales-fr@mn-net.com

USA
MACHEREY-NAGEL Inc

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

e-mail: sales-us@mn-net.com

