Certificado de calidad Quality certification Bionova® BT60/

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Steam, Ethylene Oxide and Dry Heat sterilization

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Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-2:2017, 11138-3:2017, 11138-4:2017 and IRAM 37102: 1999 (Parts 1, 2 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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Composition

Cada Bionova® BT60 envelope consists of a paper strip inoculated with a *Geobacillus stearothermophilus* ATCC 7953 and *Bacillus atrophaeus* ATCC 9372 spore population.

Product description

The envelope containing Bionova® BT60 Spore Strip is specifically designed to control Steam (Steam), Ethylene Oxide (EO) and Dry Heat (DRY) sterilization processes. Spore strips should be incubated in an appropriate growth medium for the corresponding *Geobacillus stearothermophilus* ATCC 7953 or *Bacillus atrophaeus* ATCC 9372 growth. It is suggested the use of Bionova® MC20 Culture Medium for growth of *Geobacillus stearothermophilus* ATCC 7953 at 60 °C, MC1030 for growth of *Bacillus atrophaeus* ATCC 9372 at 37 °C and/or MC1020 for growth of both spores, selective according to the incubation temperature (*Geobacillus stearothermophilus* 60 °C / *Bacillus atrophaeus* 37 °C).

If the sterilization process was not successful, MC20, MC1030 or MC1020 Growth Medium will turn to yellow after incubation at the corresponding temperature, thus indicating the existence of living spores on the strips. If the sterilization process was successful, MC20, MC1030 or MC1020 Indicating Medium will remain its original color. A final readout should therefore be performed after 24 hours of growth medium incubation at 60±2 °C (MC20 and/or MC1020) for steam sterilization processes control, or 48 hours at 37±2 °C (MC1030 and/or MC1020) for Ethylene Oxide and Dry Heat sterilization processes control.

Warning!

Do not use Bionova® BT60 Spore Strip envelopes to control radiation cycles or other sterilization processes different from those for which they have been designed. Do not reuse spore strip envelopes.

Storage

Best stored in original package under the next conditions: temperatures between 10-30 °C, 30-80 % Relative Humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Instructions for use

1. Identify the Bionova® BT60 envelope by writing the sterilizer number (in case of having more than one), load number and processing date on its back.

2. Pack it in an appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the envelope in those areas you consider a priori more inaccessible to the sterilizing agent. A typical problematic area is the center of the load.

3. Sterilize as usual.

4. After the sterilization process has finished, you can:

A. Remove the Bionova® BT60 envelope from the package for processing and incubation.

B. For Ethylene Oxide sterilization, you can aerate the package containing the biological indicator with the rest of the load or remove the biological indicator from the sterilization package for incubation prior to material aeration.

5. Tear up the envelope at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer spore strip with a sterile clamp to MC20, MC1030 or MC1020 Growth Medium tube, or to any other appropriate growth medium. **IMPORTANT:** Use latex gloves and chin strap when transferring spore strips from the envelope to growth medium tube. **WARNING!** Avoid spore strip contact with any kind of surface including the outer part of Growth Medium tube and of Bionova® BT60 envelope.

IMPORTANT: Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

6. Incubate processed spore strips along with the strip used as positive control at the corresponding temperatures: 60±2 °C or 37±2 °C. It is advisable to make observations every 10 hours.

Color change of MC20, MC1030 or MC1020 Growth Medium to yellow means a failure in the sterilization process. If after incubation, no color change is visible on the growth medium, the result is negative (i.e. the sterilization process was successful). MC20, MC1030 or MC1020 Indicating Medium used for incubating the spore strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid.

Record the positive results and discard them immediately as indicated below.

WARNING! Do not reuse the sterilizer until spore strip growth result is negative (MC20, MC1030 or MC1020 Growth Medium containing processed spore strip remains its original color).

Disposal

Discard Growth Medium tubes already used to analyse Bionova® BT60 Spore Strips growth according to your country's healthcare and safety regulations. Positive MC20, MC1030 or MC1020 Growth Medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.