

ENGLISH VERSION CATALOGUE



AGENCINOX
LA SANTÉ EXIGE LA QUALITÉ


POLYCERT
ISO 9001 : 2008
ISO 13485 : 2012

cofrac
CERTIFICATION
D'ENTREPRISES
& DE PERSONNELS
ACCREDITATION
IN ACCORD
WITH
ISO 9001 & 13485
RECOGNIZED BY
COFRAC

DC689AL



CONCEPTION, DEVELOPPEMENT, INDUSTRIALISATION

Agencinox, with a long experience in the medical field, from the medical ware to the medical furniture, offers a complete range of high quality products esteemed by the medical profession.

MEDICAL FURNITURE

Designed while bearing in mind both quality and budget aspects, Agencinox offers a range of medical furniture for the health care market: hospitals, clinics, general practitioners, laboratories...

This range of medical furniture is made of reliability and strength. Examination couch, dressing trolleys, infusion stand...

AGENCINOX : THE GUARANTEE OF QUALITY

*All items are produced in **18/10 stainless steel AISI 304** which gives the best protection against corrosion and rust thanks to its nickel content. In addition to this thorough selection of stainless steel, a high grade mirror polish finishing has been chosen for a perfect cleaning and disinfection. The combinaison of this two criterious gives a better reaction to aqueous solutions, standard detergents and diluted organic acids in the medical environment.*

CLEANING

The washing and cleaning can be done manually or in washer.

A wide range of cleaning products and methods can be used :

- *detergents, washing powders and soaps,*
- *non abrasive powders*
- *solvents*

However, we recommend the following : washing in water with a dishwashing product and rinsing.

DO NOT USE

- *Chlorinated or iodized products end their by-products*
- *Abrasive products*

The disinfection and/or sterilisation can be done by any methods such as chemical disinfection, dry heat, steam...



STAINLESS STEEL TROLLEYS

GENERAL FEATURES

Totally made of 18/10 stainless steel, bright polish finishing. Heavy duty structure. Easy cleaning.

With 2 handles, shaped trays.

Frame :

18/10 stainless steel tube diam. 25 mm, castors diam. 75 mm.

2 MODELS :

Knockdown construction or argon welded.

KNOCKDOWN CONSTRUCTION 2 OR 3 TRAYS TROLLEYS

Shaped trays. Height from floor to trays :
2 trays : 300 and 820 mm. - 3 trays : 200, 510 and 820 mm
Bumpers diam. 95 mm.

2 TRAYS

WITHOUT GUARD RAILS		WITH GUARD RAILS	
REF.	SIZE	REF.	SIZE
441	600 X 430 mm	445	600 X 430 mm
741	750 X 500 mm	745	750 X 500 mm
941	900 X 600 mm	945	900 X 600 mm

3 TRAYS

WITHOUT GUARD RAILS		WITH GUARD RAILS	
REF.	SIZE	REF.	SIZE
446	600 X 430 mm	447	600 X 430 mm
746	750 X 500 mm	747	750 X 500 mm
946	900 X 600 mm	947	900 X 600 mm

3 TRAYS DRESSING TROLLEYS

EQUIPED WITH GUARD RAILS AND STAINLESS STEEL WASTE BIN 12 L

3 removable and height adjustable rails.

Removable waste bin holder, left and right directionable. Overall size 1080x570 mm

REF.	SIZE	CASTORS
752	750 X 500 mm	100 mm
752S	Argon welded version	

2 TRAYS TROLLEYS ARGON WELDED

WITHOUT GUARD RAILS		WITH GUARD RAILS	
REF.	SIZE	REF.	SIZE
441S	600 X 430 mm	445S	600 X 430 mm
741S	750 X 500 mm	745S	750 X 500 mm
941S	900 X 600 mm	945S	900 X 600 mm

3 TRAYS TROLLEYS ARGON WELDED

WITHOUT GUARD RAILS		WITH GUARD RAILS	
REF.	SIZE	REF.	SIZE
446S	600 X 430 mm	447S	600 X 430 mm
746S	750 X 500 mm	747S	750 X 500 mm
946S	900 X 600 mm	947S	900 X 600 mm

GUARD RAILS

REF.	SIZE
441 E6	for length 600 mm
741 E6	for length 750 mm
941 E6	for length 900 mm
441 E5	for width 400 mm
741 E5	for width 500 mm
941 E5	for width 600 mm
P27	clip for guard rail

TROLLEYS

TRAYS TROLLEYS ARGON WELDED TURNED DOWN RIM SHELVES

Made of 18/10 stainless steel, argon welded.
Castor Ø100mm. Delivered without accessories.

2 TRAYS		3 TRAYS	
REF.	SIZE	REF.	SIZE
451T	600 X 400 mm	456T	600 X 400 mm
751T	700 X 480 mm	756T	700 X 480 mm



TRAYS TROLLEYS ARGON WELDED RAISED RIM SHELVES

Made of 18/10 stainless steel, argon welded
Castor Ø100mm. Delivered without accessories.

2 TRAYS		3 TRAYS	
REF.	SIZE	REF.	SIZE
451R	600 X 400 mm	456R	600 X 400 mm
751R	700 X 480 mm	756R	700 X 480 mm



Raised rim shelf



Turned down rim

TRAYS TROLLEYS ARGON WELDED TURNED DOWN RIM SHELVES

Made of 18/10 stainless steel, argon welded.
Castor Ø100mm. With guards rails.

2 TRAYS		3 TRAYS	
REF.	SIZE	REF.	SIZE
435SL	600 X 400 mm	437SL	600 X 400 mm
735SL	700 X 500 mm	737SL	700 X 500 mm
835SL	800 X 580 mm	837SL	800 X 580 mm



Ref. 737SL

TRAYS TROLLEYS ARGON WELDED RAISED RIM SHELVES

Made of 18/10 stainless steel, argon welded.
Castor Ø80 mm. Delivered without accessories.

2 TRAYS	
REF.	SIZE
351	400 X 325 mm
3 TRAYS	
REF.	SIZE
356	400 X 325 mm

OPTION : Serum rod
with clamp

REF.	DESIGNATION
609I	2 stainless steel hooks

REF.	DESIGNATION
609P	2 Polyamide hooks



Ref. 351

OPTIONS

Stainless steel drawers

FOR RAISED RIM SHELVES		TURN DOWN RIM SHELVES	
REF.	SIZE	REF.	SIZE
163	L240 X P370 X H85 mm	161	L240 X P370 X H85 mm

Stainless steel end-shelf

FOR TROLLEY		FOR TROLLEY	
REF.	SIZE	REF.	SIZE
4471	380 X 325 mm	7571	450 X 325 mm

RF	2 castors with brake
RAS	4 antistatic castors Ø 100 mm
PC	4 ring-shaped bumpers
112	Holder and basin Ø 200mm
113	Holder with tray
114	Swivel plastic bag holder 10 L
120	Swivel plastic bag holder 30 L
123	Container holder



Ref. 431N with RAS et PC



Ref. 436B with RAS et PC

TRAYS TROLLEYS ARGON WELDED WITH "COMPACT" SHELVES

Made of 18/10 stainless steel, argon welded with "compact labo tache-clean" shelves.
Castor Ø100mm. Delivered without accessories.

2 TRAYS WITHOUT GUARD RAILS		3 TRAYS WITHOUT GUARD RAILS	
REF.	SIZE	REF.	SIZE
431	600 X 400 mm	436	600 X 400 mm
731	700 X 480 mm	736	700 X 480 mm



Aigue Marine



Bleu Oslo



Nectarine



Pamplemousse

Color option, indicate after the trolley reference the color requested :

A for Aigue marine (green)

B for bleu Oslo (blue)

N for Nectarine (orange)

P for Pamplemousse (yellow)

OPTIONS

Stainless steel end-shelf

FOR TROLLEY 600x400

FOR TROLLEY 700x48

REF.	SIZE	REF.	SIZE
4471	L.380 X larg.325 mm	7571	L.450 X larg.325 mm

RF 2 castors with brake

RAS 4 antistatic castors Ø 100 mm

PC 4 ring-shaped bumpers

112 Holder and basin Ø 200mm

113 Holder with tray

114 Swivel plastic bag holder 10 L

120 Swivel plastic bag holder 30 L

123 Container holder

Stainless steel drawer

FOR TROLLEY WITH COMPACT TRAY

REF.	SIZE
165	240 X 370 X 85 mm



Ref. 113



Ref. 123 (without box)

OPTIONS

FOR ALL TYPES OF TROLLEYS



■ Holder and stainless steel waste bin Réf.110



■ Gas holder Ref.111



■ Twin bags holder 10 L Ref.220



■ Holder with basin
Ø 200 mm 2 liters
Ref.112



■ Holder with tray
260x180 mm
Ref.113



■ Swivel plastic bag holder 30 L Ref. 120



■ Guard rail and clip
ref. P27
(See page 3)



■ Swivel plastic bag holder Ref.114



■ Stainless steel drawers for shaped trays
Easy to set on 600x400 or 750 500 trolleys
Two small models can be set side by side.

Small model usefull dim.255x400x115 mm
Ref.151 (for 600 x 400) Ref.153 (for 750 x 500)

Large model usefull dim.550x400x115 mm
Ref.152 (for 600 x 400) Ref.154 (for 750 x 500)



■ Drop bottles rack
For trolleys 750x500mm
Ref. 752P



■ Probe/catheter holder
Ref.508



591CPC

ADJUSTABLE STOOLS

Seat Ø330, 60mm thick, high density foam 40kg/m³, flame retardant and washable upholstery (M1 class). Foot operated height adjustment by gas spring from 600 to 700mm. Five legs base Ø535mm, made in chromed steel, 5 castors Ø50mm with non marking type. Round tray Ø350mm for height adjustment by foot.

Seat color : ref. 591CP : Black – ref. 591CPMA : Mandarine – ref. 591CPMI : Mirabelle – ref. 591CPC : Blue Cobalt – ref. 591CPE : Blue Egee.

ADJUSTABLE STOOLS

Seat vinyl without seam Ø 330 mm. Anti-fire quality M2. Gas spring handfull adjustment (5 years warranted covered). 5 twin-castors 50 mm diameter or plastic skids.

5 legs aluminium base

REF.	Seat	Height	Base
592	330 mm	520 à 700 mm	520 mm
592CP	foot control ht. adjustable		
592D	with backrest		
592DCP	with backrest, foot control		
592P	5 plastic skids		
592PCP	5 plastic skids, foot control ht adjustable		
592DP	with backrest + 5 plastics skids		
592DPCP	backrest + 5 plastics skids, foot control ht adjustable		

5 legs nylon base

REF.	Seat	Height	Base
593	330 mm	520 à 700 mm	520 mm
593CP	foot control ht. adjustable		
593D	with backrest		
593DCP	with backrest		
593P	5 plastic skids		
593PCP	5 plastic skids, foot control ht adjustable		
593DP	with backrest + 5 plastics skids		
593DPCP	with backrest + 5 plastics skids, foot control ht adjustable		

Adjustable foot rest for stool (not for foot adjustment)

595 adjustable in height

5 legs stainless steel base. Hydraulic foot adjustment

REF.	Diam.Seat	Height	Base
598CP	330 mm skids	500 à 700 mm	520 mm
598DCP	with backrest		
598CPR	330 mm castors	520 à 720 mm	520 mm
598DCPR	with backrest		



592



592D+595



593



592DCP



598CP



REVOLVING SADDLE CHAIR FOOT OPERATED

Saddle design and comfortable seat with pommel. 40mm thick, expanded vinyl cover, washable and flame retardant (M2). Black coloured. Reclining seat to the front by socket joint. Adjustable height by gas spring from 640 to 780mm. Five legs base, Ø560mm, bended steel tube construction without any sharp edge, chrome plated. Round tray Ø350mm for height adjustment by foot. Saddle mounted on 5 rubber skids Ø45mm.

REF	DESCRIPTION
599OSC	Expanded vinyl seat, gas spring, chrome plated base on rubber skids.

SADDLE ADJUSTABLE SEAT

REFERENCE	
596	Mounted on Nylon base with 5 twin castors
596CP	With foot height adjustable
596P	Mounted on Nylon base with 5 plastics skids
596PCP	Mounted on Nylon base with 5 plastics skids, With foot height adjustable
597	Mounted on aluminium base with 5 twin castors
597CP	With foot height adjustable
597P	Mounted on aluminium base with 5 plastics skids
597PCP	Mounted on aluminium base with 5 plastics skids, With foot height adjustable



Ref. 599

SURGEON STOOL

Same as saddle but equipped with an ADJUSTABLE SITTING ANGLE allowing a SIT-STAND POSITION.

Base totally made in 18/10 aisi 304 stainless steel with foot control height adjustment from 670 to 870, this model is particularly intends for surgical operation. Mounted on rubber skids

REF	DESCRIPTION
599	Black simulated leather, hydraulic foot adjusment Stainless steel base and rubber skids.



Ref. 596

Ref. 597+595

SADDLE ADJUSTABLE SEAT

Ergonomic seat extremely comfortable made for examination and work situation. Particulary intends for generalist, specialist, dentist, ophtalmologist, gyneacologist, physiotherapist, laboratorist...

This revolutionary seat drastically reduce work fatigue due to a better blood draining.

Black simulated leather. Gaz spring handfull adjustment 5 years warranted.

Mounted on Nylon or aluminium base, with 5 twin castors or skids



OVER-BED TABLE

18/10 stainless steel. Adjustable from 890 to 1300 mm.
Authorized load 20kg.
4 antistatic castors 50 mm diameter.
Delivered knocked-down

REF.	Length	width
780	900 mm	600 mm



OVER-BED TABLE

18/10 stainless steel. Adjustable from 890 to 1300 mm.
4 antistatic castors 50 mm diameter. Authorized load 15kg
Delivered knocked-down

REF.	Length	width
781	900 mm	600 mm



MAYO TABLE

18/10 stainless steel. Adjustable from 830 to 1300 mm.
12 kg authorized load. Delivered knocked-down
4 antistatic castors 50 mm diameter.

REF.	Length	width
758	600 mm	430 mm
758F	600 mm	430 mm 2 brake wheels
758T	600 mm	430 mm Turned down shelf



HYDRAULIC MAYO TABLE

18/10 stainless steel. Adjustable from 830 to 1300 mm.
20 kg authorized load. Delivered knocked-down
4 antistatic castors 80 mm diameter.

REF.	Length	width
758H	600 mm	430 mm Embossed shelf
758LH	600 mm	430 mm Smooth shelf
758GPH	750 mm	500 mm Embossed shelf
758GPLH	750 mm	500 mm Smooth shelf



WARD SCREEN

Washable screen, classification anti-fire quality M1. Each fold is attachable to each other.


Castors diam. 50 mm. Dim. 1780x590 mm

- White screen made in polyester PVC washable, anti-fire quality M1:

Ref. 3503

- Color screen PVC, anti-fire quality M2:

Ref. 3503CH : Champagne, **3503GR** : Grège, **3503J** : light yellow, **3503A** : Amande **3503T** : Turquoise, **3503C** : Caraïbe, **3503G** : light grey



WALL WARD SCREEN

Telescopic and rotating. Adjustable length from 1120 to 2000 mm.

Washable screen, classification anti-fire quality M1. Open out length 2000 x height 1700 mm.

- White screen made in polyester PVC, washable, anti-fire quality M1

Ref. 3403

- color screen Trévira, anti-fire quality M1:

Ref. 3403S : Safran, **3403B** : blé, **3403M** : Mandarine, **3403P** : Pervenche, **3403V** : green.



Option color screen PVC anti-fire quality M2



WALL WARD SCREEN FOLDABLE

Telescopic, foldable and rotating. Adjustable length from 1200 to 2000 mm.

Washable screen, classification anti-fire quality M1. Open out length 2000 x height 1700 mm

- White screen made in polyester PVC, washable, anti-fire quality M1


Ref. 3404

- color screen Trévira, anti-fire quality M1:

Ref. 3404S : Safran, **3404B** : blé, **3404M** : Mandarine, **3404P** : Pervenche, **3404V** : green.



Option color screen PVC anti-fire quality M2



EXAMINATION LIGHT

REF. 574 20 W LIGHT

Examination halogene light 20 W on castors. Power 23000 lux/0,50m. 230/12V / 20W with protective glass. Dichroïque source (2/3 less heat). Gosse neck arm 650 mm, Head lamp : mounted on a ball joint. 2 strips for electric cable (3 m length). Nylon base with 5 castors Ø 50 mm, 18/10 stainless steel tube.

REF. 575 VARIABLE LIGHT INTENSITY 50W

Examination halogene light 50 W on castors. Power 50000 lux/0,50m. Variable intensity with dimmer switch. 230/12 V / 50 W with protective glass. Dichroïque source (2/3 less heat). Head lamp : mounted on a ball joint. Gosse neck arm 650 mm, 2 strips for electric cable (3 m length). Nylon base with 5 castors Ø 50 mm, stainless steel tube



Ref. 574



Ref. 575

REF. 576 MOBILE MAGIFYING LAMP

Electric magnifying lamp 230 V/ 22 W CE on castors. 2 articulated arms (2x456 mm). Protection IP20, magnify by 2. Two strips for electric cable setting (1m70 L). Nylon base weighted with 5 castors Ø 50 mm, 18/10 stainless steel tube.

REF. 577 MAGNIFYING LAMP

Electric magnifying lamp 230 V/ 22 W. 2 articulated arms (2x456 mm), protection IP20, magnifying by 2. Delivered with a clamp for working surface setting.



Ref. 577



Ref. 576



STEPS

Totally made of 18/10 Stainless steel.
Anti skidding steps.

REF.	MODELE	DEPTH	WIDTH	HT
759	1 step	300 mm	340 mm	210 mm
769	2 steps			
	1 st step	430 mm	340 mm	180 mm
	2 nd step	235 mm	340 mm	370 mm
769D	2 steps knocked-down construction			
779	3 steps			
	1 st step	430 mm	340 mm	180 mm
	2 nd step	333 mm	340 mm	370 mm
	3 rd step	235 mm	340 mm	560 mm
779D	3 steps knocked-down construction			



STEP ONE THREAD

Totally made of 18/10 Stainless steel.
Anti skidding step.

REF.	LENGHT	WIDTH	HEIGHT
749	600 mm	400 mm	150 mm

Special size request.

WASTE BIN

REF.	Capacity	Diam.	HT
3015	14 liters	250 mm	400mm

Stainless steel aisi 304 body, and plastic inside bin.

TOUCHLESS WASTE BIN

Top opening touchless waste bin, operated by magnet sensor which activate the lid opening. (To be used with non included battery) On/Off switch button. Stainless steel body, Removable plastic bucket with handle.

REF.	Capacity	Height	Thickness
3017	60 liters	680 mm	355mm




Ref. 3015



HEIGHT GAUGE

REF.	HEIGHT
------	--------

701 Adult height gauge 2 meterS



STRETCHER CARRIER

STRETCHER CARRIER

Totally made of 18/10 S.S. Tubes 30 mm diam.

Heavy duty stretcher ; Easily knocked down for storage. Large castors for easy moving (diam. 125 mm) with brake. Simulated leather upholstery 40 mm thickness, easy cleaning. X-rays transparent. Delivered knocked-down



986F	Overall sizes	Useful sizes
------	---------------	--------------

Length	1940 mm	1800 mm
Width	600 mm	520 mm
height	850 mm	



OPTIONS

Accessories for stretcher are adaptable at any time.

Movable side rails

REF	CARACTERISTIC
986R	Right and left, quick locking and unlocking With one hand.

Serum holder

REF	CARACTERISTIC	DIAM	length
986S	Rod, Stainless Steel 18/10 With 2 polyamide hooks and aluminium clamp.	16 mm	1000 mm

Gas holder

REF	CARACTERISTIC	size
785	18/10 stainless steel for small bottle oxygen.	Diam. 150 mm Usefull height 240 mm



EXAMINATION COUCH

REF. 585

EXAMINATION COUCH

Totally made of 18/10 Stainless steel.

Tubes Ø 30 mm.

Heavy duty stretcher. Simulated leather upholstery 40 mm thickness, X-rays transparent mattress without seam. Equipped with paper roll holder (without paper roll).

Automatic head rest locking adjustment system.

Delivered knocked-down.

REF.	Length	Width	Height
585	1930 mm	600 mm	800 mm



INSTRUMENTS TABLE 2 SHELVES ROUNDED CORNERS

Totally made from stainless steel 18/10 aisi 304 grade. Square tube frame 25 x 25 mm argon welded. 100mm antistatic castors, 4 bumpers, rounded corners, fold under shelf is 45° to facilitate the cleaning and evacuation of liquids.

Height from the floor : Lower shelf is 335mm and top shelf is 815mm.

Max load is 80kg on the 2 shelves.

Knocked down construction and 1 top shelf is available

REF.	REF.	LENGTH	WIDTH
Welded	Knocked-down		
661	663	600 mm	400 mm
761	763	750 mm	500 mm
1061	1063	1000 mm	600 mm
1261	1263	1200 mm	600 mm

INSTRUMENTS TABLE 2 SHELVES SQUARE CORNERS

Totally made from stainless steel 18/10 aisi 304 grade. Square tube frame 25 x 25 mm argon welded. 100mm castors, square corners, fold under shelf is 45° to facilitate the cleaning and evacuation of liquids.

Height from the floor : Lower shelf is 335mm and top shelf is 815mm.

Max load is 80kg on the 2 shelves.

Knocked down construction and 1 top shelf is available

Optional : bumpers, please specify the code **PC**, antistatic castors, please specify the code **RAS**.

For example : 664PCRAS means table with bumpers and antistatic castors.

REF	Length	Width
Welded		
664	600mm	400mm
764	750mm	500mm
1064	1000mm	600mm
1264	1200mm	600mm

INSTRUMENT TABLE 2 COMPACT COLOR SHELVES

Made from stainless steel 18/10 aisi 304 grade. Square tube frame 25 x 25 mm argon welded. 100mm castors, square corners, 2 shelves in "compact" resin, 4 colors. Height from the floor : Lower shelf is 335mm and top shelf is 815mm.


Max load is 80kg on the 2 shelves.

Knocked down construction and 1 top shelf is available

Optional : bumpers, please specify the code **PC**, antistatic castors, please specify the code **RAS**.

For example : 665APCRAS means table with bumpers and antistatic castors.

REF				Leng.	Width.
Aigue marine	Bleu oslo	Nectarine	Pamplemousse		
665A	665B	665N	665P	600x400mm	
765A	765B	765N	765P	750x500mm	
1065A	1065B	1065N	1065P	1000x600mm	
1265A	1265B	1265N	1265P	1200x600mm	






BABY COT

REF. 584 :

Stainless steel structure, with 4 wheels Ø50mm 2 with brakes, acrylic basket. Matress zip 50 mm washable. Delivred knocked down.. Dim. : 900 x 430 x 980 mm



WASH HAND BASIN STAINLESS STEEL

WASH HAND BASIN GENERAL FEATURES

Totally made of 18/10 stainless steel.

Leg or electronic control water.

Delivered with mixer tap.

With or without back board.



REF. _____

XB14 Leg operated, without back board

XB15 Leg operated, + back board + soap dispenser

XB35 Electronic control + back board + soap dispenser

LIZA E Round wash hand basin , Electronic control Lithium battery 6 V



**We offer a wider range of stainless steel
WASH HAND BASIN in the catalog**



ASK FOR IT!

WASH-HAND BASINS

Medical Trolleys



MEDICAL TROLLEYS

The medical trolleys are designed to meet the exact requirements of the carrying staff. They are most intended for patients care, visits patients records...

Their main characteristics are :

- Toughness (50 kg load allowed)
- Ease of cleaning
- Versality
- Attractiveness


3 SIZES

400 x 400 mm
600 x 400 mm
750 x 500 mm

Height of trolley 980 mm, top shelve 840 mm



18/10 Stainless steel shelve



2 VERSIONS OF SHELVES



2 shelves :
Height 200 mm and 840 mm

3 shelves :
Height 200, 520 and 840 mm

2 shelves or 3 shelves

MEDICAL TROLLEYS


4 OPTIONAL COLOURS



Réf. 4311M Rouge presented with accessories



Réf. 4211M Bleu presented with accessories



4 TYPES OF CASTORS (noiseless and non-marking) :



Ø 100 x 4



Ø 100 x 4
2 with brake



Ø 100 x 4
antistatic



Ø 100 antistatic
2 with brake

Packing : Delivered assembled or knocked-down in individual cardboard

☞ Possibility of specific adaptations : a large range of accessories can be fitted on these trolleys (guard rails, waste-bin, drawers...) For more details, refer to pages 23, 24.



Nb SHELVES	TYPE SHELVE	Castors Ø 100	400 x 400	600 x 400		750 x 500	
			Ref. Knocked-down	Ref. Assembled	Ref. Knocked-down	Ref. Assembled	Ref. Knocked-down
 2 Shelves	18/10 Stainless steel	4.....	3111D	4211M	4211D	5211M	5211D
		4 (2 With brake)	3112D	4212M	4212D	5212M	5212D
		4 Antistatic.....					
		4 Antistatic (2 With brake)	3113D	4213M	4213D	5213M	5213D
			3114D	4214M	4214D	5214M	5214D
 3 Shelves	18/10 Stainless steel	4.....		4311M	4311D	5311M	5311D
		4 (2 With brake)		4312M	4312D	5312M	5312D
		4 Antistatic.....					
		4 Antistatic (2 With brake)		4313M	4313D	5313M	5313D
				4314M	4314D	5314M	5314D

EMERGENCY TROLLEY 700 X 445



■ Ref. EV427

Trolley fully equipped with :

- 4 antistatic Ø 125 castors, 2 with brakes
- A work surface on upper part (at a height of 1005 mm).
- 1 heart board
- A series of drawers made up of :
 - 2 stainless steel drawers 75 mm
 - 2 stainless steel drawers 150 mm
 - 1 stainless steel drawer 230 mm

Model presented with serum rod

ANAESTHESIOLOGY TROLLEYS 700 X 445



■ Ref. EV527

Trolley without upper shelf fully equipped with :

- 4 antistatic Ø 125 castors, 2 with brakes
- A work surface on upper part (at a height of 1005 mm).
- A series of drawers made up of :
 - 1 stainless steel drawer 75 mm
 - 1 stainless steel drawer 150 mm
 - 2 stainless steel drawers 230 mm

■ Ref. EVP528

Fully equipped trolley with upper-shelf

Basically identical to ref. EV 527 but with :

- A medium + top shelf for medical apparatus

Model ref. EVP528 presented with accessories

■ VIDEO TROLLEY 5 SHELVES

The video trolley have been designed to answer to the needs of new medical technologies using video : Positioning of processor, printer, recorder, electronic bistoury, storage areas...

Dimensions : 690 x 445 mm, height 1655 mm

- 5 Stainless steel shelves 620x400 mm
- 4 antistatic castors Ø 100, 2 with brake

Ref. 4514M (assembled)

Ref. 4514D (knocked-down)



EMERGENCY TROLLEYS

MEDICAL TROLLEYS WITH STORAGE UNIT

THE EVOLUTION RANGE

This range has been designed to match the specific requirements of each hospital (I.C.U., emergency...)

Two versions :

Ref. **EV** : Basic trolley (without drawers)

Ref. **EVP** : Basic trolley with a gantry (without drawers)



Model : EV + 1B/1B/2B/2B/3B

CHARACTERISTICS :

- Work surface made in 18/10 polished stainless steel HT 1005 mm
- Overall size 700 x 445 x 1140 mm
- 4 antistatics wheels Ø 125 (2 with brake)+ 4 plastics bumper

➤ Sturdiness :

- Square tube 25 x 25 x 1,2 polished stainless steel 18/10
- Insert for wheel welded inside the tube.
- Carrying capacity : 100 kg

➤ Easy to clean : top surface without gap.

➤ Easy to move : 4 antistatics wheels 125 mm Ø

➤ Secure : central locking cross bar (optional)

➤ Ergonomic : large drawer (600 mm) full opening.

➤ Modular : multi purpose, personalized, evolution.

Trolleys (EV or EVP) are designed to receive 9 units of drawers

3 differents heights :

- | | |
|-----------------------------|----------------------------------|
| ● 1 drawer 75 mm = ① unit | 1 drawer 75 mm height : 1 unit |
| ● 1 drawer 150 mm = ② units | 1 drawer 150 mm height : 2 units |
| ● 1 drawer 230 mm = ③ units | 1 drawer 230 mm height : 3 units |

Drawers available in two versions :

- Stainless steel front
- coloured "Compact" front :



Aigue Marine



Bleu Oslo



Nectarine



Pamplemousse

Model EVP + 1B/2N/3B/3N Options : drawers locking
EV C, EVPB5, tray : EVP T and EVP E



"EVOLUTION"
MEDICAL TROLLEYS

● Front panel references :


- Stainless steel drawer front panel: **réf. I**
- Nectarine drawer front panel : **réf. N**
- Blue drawer front panel : **réf. B**
- Yellow drawer front panel : **réf. P**
- Green drawer front panel : **réf. A**

COMPOSITIONS FOR EXAMPLE


For creating your trolley with storage unit select the trolley frame (EV or EVP)
then add the drawers

All stainless steel and coloured Compact front panel drawers can be used on the same trolley.

Important : When adding the drawers you have to obtain 9 units



Trolley **EV** + 1N / 1N / 2B / 2B / 3N
For the gantry, precede the ref. **EVP**



Trolley **EV** + 1P / 2B / 2B / 2B / 2P

OPTIONS : Options are the same as for medical trolleys :

- **ref. EVC** : Central locking cross bar
- **ref. 650** : Heart board 600x390 ép.10mm.
- **ref. EVG** : Guard rail for EV and EVP.
- **ref. EV742** : Swivel monitoring tray 430x360 mm.
- **ref. EV71** : Stainless steel retractable side shelf.
- **ref. 707** : 5 Multi sockets panels

FOR EVP

Basic model with a gantry.

- **ref. EVP E** : Top shelf
- **ref. EVP T** : Medium shelf
- **ref. EVP B6** : Block with 6 swinging plastic boxes and stainless steel holder
- **ref. EVP B5** : Block with swinging plastic boxes and stainless steel holder
- **ref. EVP B4** : Block with 4 swinging plastic boxes and stainless steel holder



Ref. EV742

DRESSING AND TRANSPORTING TROLLEYS



Ref. 655 Model presented with accessories

Stainless steel removable rails
Possibility to set each rail every 75 mm
Castors Ø 100 mm, 2 with brake + bumpers.

Rails x 2 for trolley ref. 655

Ref. 6604

SINGLE TROLLEY Ref. 655

Length 655 x width 450, height 1020 mm
Delivered without rails and trays.

STAINLESS STEEL WORKING SURFACE Ref. 655P

For single trolley

ACCESSOIRES FOR ALL CARTS

CATHETER HOLDER
Adaptable on medication cart

Ref. 509

LARGE BAG HOLDER 30 L

Ref. 620

TRAY HOLDER
with tray 260 x 180 x 30 mm

Ref. 604

TWIN BAGS HOLDER 10L

Ref. 720

BASIN HOLDER
with basin 2 liters Ø 200 mm

Ref. 605


BAG HOLDER 8L

Ref. 606

MULTI SOCKETS PANEL 5 CONNECTIONS

Ref. 607

POLYCARBONATE TRAYS AND DIVIDERS



International standardized sizes : 600 x 400 mm
3 heights : 50, 100, 200 mm
Tinted polycarbonate trays, can be washed in high temperature machine and resistant to standard detergents.
Possibility for dividers in length and width.

POLYCARBONATE TRAYS

600 x 400 x 50 mm
600 x 400 x 100 mm
600 x 400 x 200 mm


Ref. 6050

Ref. 6100

Ref. 6200

- Tray height 50 mm could be adjusted every 75 mm
- Tray height 100 mm could be adjusted every 150 mm
- Tray height 200 mm could be adjusted every 225 mm

Possibility to combine trays each others



DIVIDERS FOR TRAYS

Length divider 600 x 50 x Ht 50mm

Ref. 6052

Width divider 400 x 50 x Ht 50 mm

Ref. 6053

Length divider 600 x 100 x Ht 100mm

Ref. 6102


Width divider 400x100 or length for ref. 6100,

Ref. 6103

ACCESSORIES

Modular, easy to fit, they offer great flexibility of use and adaptation.

GUARD RAILS



For trolley 600 x 400

Ref. 460 For length
Ref. 440 For width

For trolley 750 x 500

Ref. 560 For length
Ref. 540 For width

Sold in single units with fixing. Bear in mind :
Number of guard rail per shelf.

14L WASTE BIN

Ref. 600

Stainless steel 18/10
With holder and fixing



DRAWERS



Made in stainless steel. Quickly set up. Installed on front face

Small drawer
602P for trolley 600x400
602G for trolley 750x500

Useful dimensions :
370 x 240 x 85 mm

Large drawer
603P for trolley 600x400
603G for trolley 750x500

370 x 535 x 85 mm

TRAY HOLDER

Ref. 604

Stainless steel with
Tray 260x180x30
260 x 180 x 30 mm
Pivoting holder



PLASTIC BAG HOLDER - Ref. 606




Pivoting, made in
stainless steel

Ref. 616

50 plastic bags 8 L

RAILS



Delivered with fixations

For trolley 600 x 400

Ref. 403 For length

For trolleys 750 x 500 :

Ref. 503 For length

SERUM ROD



Adjustable on Agencinox trolleys with
clamp

Ref. 609 P
Rod with 2 Polyamid hooks


Ref. 609 I
Rod 2 stainless steel hooks



BASIN HOLDER

Ref. 605

Stainless steel with basin
capacity 2 L, Ø 200 mm
Pivoting holder



MULTI SOCKETS PANEL Ref. 607 FOR 5 CONNECTIONS



Fixed under the edges
of the shelf.
Switch with light
indicator

ACCESSORIES

ACCESSORIES



Ref. 608
GAS CYLINDER HOLDER

Stainless steel with fixing
pivoting holder for medical
gas cylinder Ø 150 mm



Ref. 509
PROBE CATHETER HOLDER

Stainless steel with fixing
Dim. 160 x 70 x 400 mm



BAG HOLDER 30 L

Made in 18/10 stainless
steel

Ref. 620 :



TWIN BAG HOLDER
10 L


Made in 18/10 stainless
steel

Ref. 720



Ref. 645
DROP BOTTLE BRACKETS

For use directly on guard rails
For all diameter of bottles



RETRACTABLE END-SHELF

Ref. 471 : For trolley 600 x 400
Dim. 376 x 300 mm

Ref. 571 : For trolley 750 x 500
Dim. 456 x 300 mm

PARTITION FOR DRAWERS (BY 10)



REF	Length	Width	Height
614	150 mm	50 mm	43 mm
613	75 mm	50 mm	43 mm
615	225 mm	50 mm	43 mm
617	100 mm	75 mm	43 mm
618	150 mm	100 mm	43 mm

ISOLATION CART

NEW !



Ref. 982 ISOLATION CART

Totally made from stainless steel, square tube 25x25mm.

Easy to move, modular and evolutionary, this cart provide a quick isolation and protection before entering in the patient room.

Cart specifications :

- 1 Stainless steel square tube frame Height 1740mm, W.450mm, D.610mm.
- 1 White board screen for instructions writing.
- 1 Alcohol or gel flask holder.
- 1 Clothing holder.
- 4 swiveling antistatic castors Ø50mm including 2 with brake.

Reduced overall dimensions by the optimal size and nice design, it is customizable with the following options :

Options shown on the right pictures :

- | | |
|---|-------------------|
| ➤ Wire basket St. St. 480 x 300 x 150mm | ref. ST63 |
| ➤ St. St. wire shelf 480 x 380 x 150mm | ref. ST64 |
| ➤ St. St. sheet shelf 480 x 380 x 150mm | ref. ST65 |
| ➤ Epoxy coated wire basket 7035 480x300x150mm | ref. ST69 |
| ➤ St. St. wire basket (inclined) | ref. ST70. |
| ➤ Revolving bag holder right side | ref. 114D. |
| ➤ Revolving bag holder left side | ref. 114G |
| ➤ Container holder (without box) | ref. 124 |

NEW !

EASY CLEANING BASE

STAINLESS STEEL HOOKS





**STAINLESS STEEL
LOCKING SLEEVE**

OPTIONAL FUNNY CASTORS FOR CHILDREN



EASY CLEANING




Description

Serum infusion stand with plastic base on wheels (500mm width base) easy to clean, without harshness or cavity. Available in 6 colors, black, grey, blue, red, yellow and green. Lower tube and rod telescopic sleeve Ø25mm Ø16mm 18/10 stainless steel AISI 304. Locking of the telescopic rod with brass screw on robust steel sleeve. Maximal capacity load 2 kg per hook.

Ref.	Description	Height maxi
1773 P	With 2 Polyamid hooks	2000mm
1773 I	With 2 stainless steel hooks	2000mm
1773 4	With 4 stainless steel hooks	2000mm
1773 U	With 2 stainless steel hooks « U »	2000mm
1773 4U	With 4 stainless steel hooks « U »	2000mm

Choose your base color by using a letter as following :

B for **Blue** / G for **Grey** / R for **Red** / J pour **Yellow** / V pour **Green**




RANGE OF INFUSION STAND

■ 5 TYPES OF STANDS AVAILABLE

Cast iron stand
Covered with 18/10 stainless steel
Ø 280 mm



4 castors stainless steel stand
Square tube 35 x 35 mm
600 mm width base



5 castors stainless steel stand
Square tube 25 x 25 mm
600 mm width base



5 castors moulded Nylon stand
520 mm width base



5 castors aluminium stand
Polished moulded aluminium
520 mm width base



■ CASTORS FOR STAINLESS STEEL STANDS

2 types of castors are available



Twin black castors
Ø 50 mm (authorized load
40 kg/castor)



Antistatic castors
Ø 50 mm Subject to regulation with regard to
gas anaesthetics in ops room

- Weighting available (screwed on the bottom of stand). Obligatory for infusion stands with a maximum height of 2400 mm.

■ CASTORS FOR NYLON AND ALUMINIUM STANDS

5 types of castors are available



Twin black castors Ø 50 mm
(authorized load 40 kg/castor)



Twin black castors Ø 50 mm (authorized load 40 kg/castor)
Associated with coloured hooks and sleeves.

- Weighting available (screwed on the bottom of stand).
Obligatory for infusion stands with a maximum height of 2400 mm.

■ STANDARD SLEEVE

For all stands


18/10 stainless steel tube Ø 30 mm –
Length 1000 mm
With stainless steel locking sleeve.




■ SERUM RODS

18/10 stainless steel Ø 16 mm tube
1000 mm length allowing height adjustment to 2000 mm
1500 mm length allowing height adjustment to 2400 mm

5 types of safety hooks



2 unbreakable polyamide hook
3 kg authorized per hook



2 stainless steel hooks
authorized load 12 kg
(6 per hook)



4 stainless steel hooks
authorized load 12 kg
(3 per hook)



4 plastic hooks
authorized load 12 kg
(3 per hook)

■ ACCESSORIES

To meet the requirements of the different departments a wide range of accessories is available at the end of this brochure (instrument trays, serum holders...)



Réf. 770 I



Réf. 770 4

■ 2 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height
770 I	2000 mm
770 IL	2400 mm

■ 4 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height
770 4	2000 mm
770 4L	2400 mm



Réf. 1774 I



Réf. 1774 P

■ 2 POLYAMIDE SAFETY HOOKS

Reference	Maxi height
1774 P with twin castors	2000 mm
2774 P with antistatic castors	2000 mm

■ 2 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height
1774 I with twin castors	2000 mm
2774 I with antistatic castors	2000 mm

■ 4 STAINLESS STEEL SAFETY HOOKS

Reference :	Maxi height
1774 4 with twin castors	2000 mm
2774 4 with antistatic castors	2000 mm



Twin black castors Ø 50 mm or
Antistatic castors Ø 50 mm
Black color



Ref. 2775 I

Ref. 2775 4

Twin black castors Ø 50 mm or
Antistatic castors Ø 50 mm
Black color



Ref. 1775 P



■ 2 POLYAMIDE SAFETY HOOKS

Reference	Maxi height	Weighted
-----------	-------------	----------

With twin castors

1775 P	2000 mm	no
3775 P	2000 mm	yes
3775 PL	2400 mm	yes

With antistatic castors

4775 P	2000 mm	yes
---------------	---------	-----

■ 2 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height	Weighted
-----------	-------------	----------

With twin castors

1775 I	2000 mm	no
3775 I	2000 mm	yes
3775 IL	2400 mm	yes

With antistatic castors

2775 I	2000 mm	no
4775 I	2000 mm	yes

■ 4 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height	Weighted
-----------	-------------	----------

With twin castors

1775 4	2000 mm	no
3775 4	2000 mm	yes
3775 4L	2400 mm	yes

With antistatic castors

2775 4	2000 mm	no
4775 4	2000 mm	yes
4775 4L	2400 mm	yes

STAINLESS STEEL STAND 5 CASTORS



Twin castors Ø 50 mm
Black color



Twin castors Ø 50 mm color

1776 P



1776 PR



1776 4V



1776 I



1776 U

■ POLYAMIDE SAFETY HOOKS

Reference	Maxi height	Weighted
1776 P	2000 mm	no
1776 4P	2000 mm	no
3776 P	2000 mm	yes
3776 4P	2000 mm	yes
3776 PL	2400 mm	yes

Colours option : indicate after the stand reference the colour requested, R for red, V for green, J for yellow, B for blue.

■ 2 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height	Weighted
1776 I	2000 mm	no
3776 I	2000 mm	yes
3776 IL	2400 mm	yes

Colours option : indicate after the stand reference the colour requested, R for red, V for green, J for yellow, B for blue.


■ 2 STAINLESS STEEL OPEN HOOKS

Reference	Maxi height	Weighted
1776 U	2000 mm	no

■ 4 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height	Weighted
1776 4	2000 mm	no
3776 4	2000 mm	yes
3776 4L	2400 mm	yes

Colours option : indicate after the stand reference the colour requested, R for red, V for green, J for yellow, B for blue.



Twin castors Ø 50 mm
Black color



Twin castors Ø 50 mm color

3778 P



■ 2 POLYAMIDE SAFETY HOOKS

Reference	Maxi height	weighted
-----------	-------------	----------

1778 P	2000 mm	no
3778 P	2000 mm	yes
3778 PL	2400 mm	yes

Colours option : indicate after the stand reference the colour requested, R for red, V for green, J for yellow, B for blue

■ 2 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height	Weighted
-----------	-------------	----------

1778 I	2000 mm	no
3778 I	2000 mm	yes
3778 IL	2400 mm	yes

Colours option : indicate after the stand reference the colour requested, R for red, V for green, J for yellow, B for blue.

■ 4 STAINLESS STEEL SAFETY HOOKS

Reference	Maxi height	Weighted
-----------	-------------	----------

1778 4	2000 mm	no
3778 4	2000 mm	yes
3778 4L	2400 mm	yes


Colours option : indicate after the stand reference the colour requested, R for red, V for green, J for yellow, B for blue

ALUMINIUM STAND 5 CASTORS



Antistatic castors Ø 50 mm
Black color

Ref. 5775 4L



■ **ÉQUIPMENT :**

16 mm diameter tube with 4 safety hooks. 1620 mm allowing maximum height adjustment to 2980 mm.
30 mm diameter standard sleeve – length 1500 mm.
5 castors stainless steel stand with antistatic and 2 locking castors 50 mm diameter.

■ **CAPACITY FEATURES :**

To adapt up to 5 syringe pumps, infusion pumps...
Multi sockets panel for 5 pieces of equipment.

■ **ERGONOMICS AND BALANCE**

Resseced castors for easy walking and pushing.
Secured balance by speading center of gravity.

■ **ACCESSORIES :** (See page 38)

- Weighting 6,5 kg
- Round tray Ø 350 mm
- Square tray 315 x 275 mm
- Multi sockets panel for 5 connections, switch with light indicator, cable length 2600 mm

Réf. 770AS

Réf. 773

Réf. 772

Réf. 607 P



SERUM ROD

18/10 stainless steel – 16 mm diameter.
Adaptable on beds, wheelchairs, trolleys, wall brackets...

■ 1000 MM LENGHT RODS

2 Polyamide safety hooks : **Réf. P**

Colours option : indicate after the stand reference the colour requested, **R** for red, **V** for green, **J** for yellow, **B** for blue.

2 stainless steel safety hooks : **Réf. I**

■ 1500 MM LENGTH RODS

2 Polyamide safety hooks : **Réf. PL**

Colours option : indicate after the stand reference the colour requested, **R** for red, **V** for green, **J** for yellow, **B** for blue

2 stainless steel safety hooks : **Réf. IL**

■ CRANKED RODS

2 stainless steel safety hooks : **Réf. IC**

4 stainless steel safety hooks : **Réf. 4C**

■ 1000 MM LENGTH RODS



2 stainless steel hooks « U » : **Réf. IU**

4 stainless steel hooks « U » **Réf. 4U**

■ 1500 MM LENGTH RODS

2 stainless steel hooks « U » : **Réf. ILU**

4 stainless steel hooks « U » **Réf. 4LU**





SERUM RODS

■ UNIVERSAL MEDICAL CLAMPS

Made of anodized aluminium

Horizontal and vertical setting :

- Rail height 40 mm maximum width 8 or 10 mm (15 kg authorized load)
- Tube Ø from 12 to 40 mm (6 kg authorized load)
- Square tube from 12 to 40 mm (15 kg authorized load)
- Square tube from 12 to 30 mm (angle setting) (15 kg authorized load)
- Shelf from 0 to 40 mm (15 kg authorized load)



Ref. 7006

For serum rods setting from Ø 7 to 16 mm
Delivered without serum rod

Ref. 7007

For medical strips apparatus

CLAMPS



■ **SERUM HOLDER BASKET -
PLASTIC COATED**

- For 250 ml bottle **Réf. F250**
For 500 ml bottle **Réf. F500**
For 1000 ml bottle **Réf. F1000**



■ **TRAYS FOR INFUSION STANDS**

■ **ROUND TRAY** **Réf. 773**

18/10 stainless steel – Ø 350 mm
Adjustable height – Locking with screw
Authorized load : 15 kg



■ **SQUARE TRAY** **Réf. 772**

18/10 stainless steel
Useful dimensions : 315 x 275
Adjustable height – Locking with screw
Authorized load : 15 kg



■ **WEIGHTING** **Réf. 770 AS**

Adaptable to all infusion stand with castors.
Ø 160 mm – Height 45 mm
Weight : 6,5 kg



■ **HANDLE** **Réf. 790**

Adaptable to all infusion stand with castors
Ø 30 mm



■ **URINE BAG HOLDER** **Réf. 774**

Adaptable to all infusion stand with castors.
Ø 30 mm



■ **MULTI-SOCKETS PANEL
FOR 5 CONNECTIONS**

Switch with light indicator,
cable length 2600 mm

Réf. 607P

■ **GAS HOLDER**

For bottle Ø 150 mm

Réf. 785



■ **SINGLE BOWL STAND**

Réf. 775

Height 860 mm – single bowl Ø 280 mm
5 antistatic castors Ø 50 mm stainless steel
stand, square tube 25 x 25 mm.





■ **DOUBLE BOWL STAND**

Réf. 778

Height 860 mm – double bowl Ø 280 mm
5 antistatic castors Ø 50 mm stainless steel
stand, square tube 25 x 25 mm.

STANDARD MEDICAL RAIL AND ACCESSORIES

Rail size complies with standard NF EN ISO 19054. This medical rail is made in anodized aluminium.




4 sizes available : 1M, 2M, 3M, delivered with wall contact elements and end stop. We can also sell medical rail at the desired length.

RAIL LENGTH	Reference	
1 meter	RAIL 100	Delivered with 2 wall contact elements and 2 end stop
2 meters	RAIL 200	Delivered with 4 wall contact elements and 2 end stop
3 meters	RAIL 300	Delivered with 6 wall contact elements and 2 end stop
Length sale	RAIL : Add after the reference RAIL the desired length in cm and order the number of wall contact elements needed and end stop as per below instructions.	

Maximum length between 2 wall contact elements :

Refer to the draw below for ordering the correspondent wall contact elements.



Setting :

It is essential to select the appropriate screws (non furnished) for fixing the wall contact elements as per the type of the wall (brick wall, breeze block, plasterboard...).

RAIL ACCESSORIES

For all our accessories the correspondent clamp is included

SERUM RODS



7015 7016 7017 7018


18/10 stainless steel. Tube Ø 16 mm
Delivered with clamp.
6kg authorized load



7018

PROBE/CATHETER HOLDER

Stainless steel suspended probe/catheter holder




507 507T

Ref. 507 : Height : 395 mm / Width : 160 mm
Depth : 70 mm / 1 compartment

Ref. 507T : Height : 510 mm / Width : 380 mm
Depth : 70 mm / 3 compartments

SWIVEL ARM AND SETTING COLUMN



Ref. 7200 :
Clamp for swivel arm or.
Column. 15 kg authorized
load.



Ref. 7201 :
Swivel arm 200 mm
10 kg authorized load




Ref. 7100 :
Column 1000 mm for clamp **7200** or
Swivel arm **7201**. These stainless steel
columns (Ø 30 mm) are made for setting
syringue pump or other apparatus
equipped with clamp.
15 kg authorized load




Ref. 7220 :
Swivel tray + clamp, height : 15 mm
Width : 340 mm, depth : 300 mm
This tray could be coupled with
a swivel arm ref. **7201**

■ **SHELVES**



Ref. 7075 :

Fix stainless steel shelf with rim. Width : 365 mm, depth : 290 mm
Carrying capacity : 15 kg. Delivered with 2 clamps.



Ref. 7575 :

Foldable stainless steel shelf without rim.
Width : 465 mm, Depth : 300 mm. Carrying capacity 15 kg.
Delivered with 2 clamps.

NOTA : Loading test capacity had been done with our rail page 37.

■ **EXAMINATION LIGHT**



Ref. 7005 :

Halogene light 12 volts / 20 watts, 9000 lux at 0,50 meter
3 articulated parts for adjustment

■ **STAINLESS STEEL CHART HOLDER**



Ref. 422 : Suspended chart holder, width 300 x height 210 mm

■ **BASIN**



Ref. 7112 : Holder with stainless steel basin Ø 200 mm (2 L) + clamp

■ **TRAY**



Ref. 7113 : Holder with stainless steel tray 260 x 180 mm + clamp

■ **MULTI-SOCKETS PANEL**



Ref. 7116 : 5 connections with light indicator + clamps

■ **SWIVEL PLASTIC BAG HOLDER**



Ref. 7117 : Holder for bag 8 L + clamp


Ref. 616 : 50 plastic bags 8 L capacity

■ **HANGING HOOK**



Ref. 7118 : Hanging hook + clamp
Possibility for hanging hook ramp or other size on request.
10 kg authorized load

Medical Ware





BED PAN WITH HANDLE (PLASTIC)

Ref.	Length	width	height	capacity
1057	450 mm	250 mm	100 mm	2,5 liters

Ref.	
1058	Lid for plastic bed pan

Ref.	
1067	Stainless steel holder for bed pan



PLASTIC URINALS

Ref.		Length	width	height	capacity
1050	for man	310 mm	140 mm	135 mm	1 litre
1053	for woman	340 mm	100 mm	120 mm	1 litre



URINE BOTTLE

Ref.	Diam.	height	Capacity
1055	125 mm	245 mm	2 litres

MEDICAL PEAR SHAPED

Ref.	Designation
412	Pear shaped 25 ml



IRRIGATOR 18/10 STAINLESS STEEL

Ref.	Diameter	height	Capacity
4023A	120 mm With handle	180 mm	1,7 L

ACCESSORIES FOR IRRIGATOR

Ref.	Designation
402	Vaginal, rectal, cock and pipe 1M50
402T	Caoutchouc pipe 1M50



TONGUE DEPRESSOR ANGLE 135° 18/10 STAINLESS STEEL

Ref.	Model	thickness	length.
101	adult	15/10 ^{ème}	160 mm
101B	baby	15/10 ^{ème}	150 mm



RECTANGULAR TRAYS FOR INSTRUMENTS - 18/10 STAINLESS STEEL

Ref.	Length	Width	Height
1217	170 mm	110 mm	20 mm
1220	205 mm	130 mm	20 mm
1222	220 mm	145 mm	28 mm
1224	245 mm	175 mm	30 mm
1226	265 mm	185 mm	30 mm
1228	285 mm	185 mm	30 mm
1230	300 mm	210 mm	40 mm
1235	355 mm	255 mm	40 mm
1242	425 mm	280 mm	50 mm
1252	510 mm	340 mm	50 mm



DENTAL TRAYS - 18/10 STAINLESS STEEL

Ref.	Length	width	Height
1210	205 mm	105 mm	15 mm
1215	205 mm	145 mm	15 mm



JARS - 18/10 STAINLESS STEEL

Ref.	Model	Diam.	Height
7011	jar thermometer	30 mm	110 mm
7030	Pean or Kocher forceps jar	30 mm	150 mm
7040	serving forceps jar	40 mm	230 mm



GALLIPOTS - LOTION BOWLS

Quantity By 5

Ref.	Diam.	Height	Capacity
210	100 mm	75 mm	400 ml
215	150 mm	75 mm	1 litre




RECTANGULAR TRAYS - AUTOCLAVABLES

Ref.	Length	Width	Height
2015	200 mm	150 mm	50 mm
3025	300 mm	250 mm	50 mm



DROP BOTTLE

Ref.	Diam.	Height	Capacity
4021	45 mm	160 mm	125 ml
4022	60 mm	200 mm	250 ml
4025	70 mm	220 mm	500 ml



COLOURED INSTRUMENTS BOXES ANODISED ALUMINIUM

Ref.	Length	Width	Height
8029 (B,R,J,V)	170 mm	70 mm	30 mm
8030 (B,R,J,V)	180 mm	90 mm	30 mm
8031	180 mm	90 mm	40 mm
8034	215 mm	100 mm	30 mm
8036 (B,R,J,V)	215 mm	100 mm	50 mm

For ordering the coloured boxes, please add the letters **B** for blue, **R** for red, **J** for yellow, **V** for green

INSTRUMENTS BOXES 18/10 STAINLESS STEEL

Ref.	Length	Width	Height
5029	170 mm	70 mm	30 mm
5030	180 mm	90 mm	30 mm
5031	180 mm	90 mm	40 mm
5036	215 mm	100 mm	50 mm
5046	250 mm	150 mm	60 mm
5052	300 mm	150 mm	70 mm
5062	340 mm	200 mm	80 mm
5069	400 mm	200 mm	90 mm
5075	450 mm	200 mm	100 mm

STAINLESS STEEL BOXES

Can be manufactured in other sizes

PERFORATED BOXES ANODIZED ALUMINIUM

Ref.	Length	Width	Height
9030	180 mm	90 mm	30 mm
9036	215 mm	100 mm	50 mm

COTTON CONTAINER 18/10 STAINLESS STEEL

Ref.	with lid	Diam.	Height	Capacity
5002		125 mm	125 mm	1 litre

Do not use any soda whashing products for all aluminium boxes.



LOTION BOWLS WITH SPOUT - 18/10 STAINLESS STEEL

Ref.	Diam.	Height	Capacity
304 by 5	40mm	20mm	0,020 L
305 by 5	50mm	25mm	0,035 L
306	60mm	30mm	0,060 L
307	70mm	35mm	0,085 L
308	80mm	40mm	0,100 L
310	100mm	50mm	0,250 L
312	120mm	60mm	0,500 L
316	160mm	80mm	1 L
320	200mm	100mm	2 L
325	250mm	120mm	4,5 L

For ordering lotion bowls without spout, please add « **SSB** »
After the reference and with handles, add « **A** » (from ref.308 to ref. 325)

BASINS AND WHASH BOWLS 18/10 S.S.

Ref	Diam.	Height	Capacity
1520	200mm	115mm	2 L
1528	280mm	125mm	5 L
4035	400mm	155mm	14 L

KICK BUCKET WITH TROLLEY 18/10 S.S. 4 CASTORS 50 MM

Ref.	trolley
4045RA	Model with antistatics castors with S.S. basin 14L

KICK BUCKET WITH TROLLEY 18/10 S.S. 3 CASTORS 50 MM

Noiseless rubber joint

Ref.	
4046RA	Model with antistatics castors with S.S. basin 14L

KIDNEY DISH 18/10 STAINLESS STEEL

Ref.	Lenght	Width	Height	Capacity
Model without rim				
1717	170mm	95mm	35mm	0,20 L
1724	250mm	135mm	45mm	0,75 L
Model halfround rim				
1725	250mm	135mm	45mm	0,75 L



BEAKERS MOLYBDENE STAINLESS STEEL

Ref.	Diam.	Height	Capacity
without spout			
2010	50 mm	60 mm	0,100 L
2025	65 mm	80 mm	0,250 L
2050	80 mm	110 mm	0.500 L
2100	100 mm	130 mm	1 litre
2200	120 mm	180 mm	2 litres
2300	150 mm	180 mm	3 litres
2500	180 mm	210 mm	5 litres

For option ordering please add after the reference **B** for spout, **A** for handle, **AB** for handle and spout.

HIGH EVAPORATING DISHES

Ref.	Diam.	Height	Capacity
1070	70 mm	35 mm	125 ml
1085	85 mm	45 mm	250 ml
1110	110 mm	55 mm	500 ml
1140	140 mm	70 mm	1 L
1160	160 mm	105 mm	2 L
1200	200 mm	160 mm	5 L

EVAPORATING DISHES

Ref.	Diam.	Height	Capacity
406	60 mm	20 mm	50 ml
408	80 mm	20 mm	100 ml

FUNNELS WITH HANDLE

Ref.	Diam.	Capacity
4215	150 mm	0,5 litre
4220	195 mm	1 litre
4250	250 mm	3 litres

SPATULAS BY 5

Ref.	Model	Length
1869	spatula	195 mm
1975	double spatula	230 mm
1822	shaft spatula	130 mm
1826	micro spatula	130 mm
1975LG	double spatula	500 mm

CHEMICAL SCOOP

Ref.	Diam.	Height	Capacity
915	75 mm	150 mm	0,5 L
920	95 mm	200 mm	1 L
923	115 mm	233 mm	2 L



OPERATING AREA AND STERILIZATION FURNITURE





AGENCINOX : THE GUARANTEE OF QUALITY

All items are produced in **18/10 stainless steel AISI 304** which gives the best protection against corrosion and rust thanks to its nickel content. In addition to this thorough selection of stainless steel, a high grade mirror polish finishing has been chosen for a perfect cleaning and disinfection. The combination of this two criterions gives a better reaction to aqueous solutions, standard detergents and diluted organic acids in the medical environment.

CLEANING

The washing and cleaning can be done manually or in a washer.

A wide range of cleaning products and methods can be used :

- detergents, washing powders and soaps,
- non abrasive powders
- solvents

However, we recommend the following : washing in water with a dishwashing product and rinsing.

DO NOT USE

- Chlorinated or iodized products and their by-products
- Abrasive products

DISINFECTION

The disinfection and/or sterilisation can be done by any methods such as chemical disinfection, dry heat, steam...

AGENCINOX

Technical benefit

Designed to adapt to your storage area, the shelving is entirely made of stainless steel pipe 18/10 square with a 25x25 mm section and equipped with support struts providing excellent *robustness*.

For better hygiene, our ladder stiles are plugged on top.

FIXED SHELVING

Equipped with adjustable pads to compensate for the unevenness of the floor, our ladders guarantee optimal stability of your shelving.

ROLLING SHELVING

Particularly easy-to-handle, the rolling shelving is equipped with wheels with a 125mm diameter including 2 with brakes for optimal safety.

In order to be cleaned in a washer, your shelving may be equipped with special wheels.

Practical benefit

In order to optimize your storage area, our ladders are available in two heights:

- | | |
|---|-----------------------------|
| > 1,700 mm consisting of 9 cross-members and as many shelves as possible | FIXED
SHELVING |
| > 2,000 mm consisting of 11 cross-members and as many shelves as possible | |
| > 1,110 mm consisting of 5 cross-members and as many shelves as possible | ROLLING
SHELVING |
| > 1,855 mm consisting of 9 cross-members and as many shelves as possible. | |

Depth :450mm or 600mm.

The cross-members of the ladders are welded at a pitch of 160 mm, leaving a sufficient working space between the shelves for all your instrument containers, usually with a height of 150 mm.

Placing one shelf per cross-member will optimize your storage.

For instrument containers higher than 150 mm, you should place a shelf every other cross-member.

In order to facilitate the layout, we suggest four shelving lengths for an even depth of 600mm or 450mm.

4 available lengths: 700 mm/1,000 mm/1,300 mm/1,500 mm

Each establishment has its own features (packing, cleaning, etc.) and, as a result, we adapt to your practices and we provide two shelving finishes:

- > stainless steel wire shelving 18/10
- > stainless steel full shelving 18/10

The modular benefit

To customize your storage area based on your premises and passage ways, we suggest modular solutions:

- > Fixed shelving used alone
- > Rolling shelving used alone
- > Fixed shelving with right angle extension or aligned

FIXED STAINLESS STEEL SHELVING

Shelving sizes	600x700	600x1000	600x1300	600x1500
Height 1700mm	RFS17070	RFS17010	RFS17013	RFS17015
Height 2000mm	RFS20070	RFS20010	RFS20013	RFS20015
Shelving sizes	450x700	450x1000	450x1300	450x1500
Height 1700mm	RFS17070.450	RFS17010.450	RFS17013.450	RFS17015.450
Height 2000mm	RFS20070.450	RFS20010.450	RFS20013.450	RFS20015.450

Consisting of : 2 ladders
 4 wire shelves for a shelving height of 1700mm
 Or 5 wire shelves for a shelving height of 2000 mm
 4 support struts



Create your own storage area

18/10 Stainless steel ladders

LADDERS	600x1700mm 9 Shelves	600x2000mm 11 Shelves	450x1700mm 9 Shelves	450x2000mm 11 Shelves
Reference	EC170	EC200	EC170.450	EC200.450

Support struts (sold in pairs)

They strengthen the structure.

- > For single fixed shelving : 2 pairs of struts
- > For extension basic shelving : 1 pair of struts
- > For aligned extension : 1 pair of struts per module
- > For herring-boat pattern : no struts.

STRUTS	700mm	1 000mm	1 300mm	1 500mm
Reference	ENT70	ENT100	ENT130	ENT150

Shelves

Designed to be directly installed on the ladder cross-members

SHELVES	600 x 700	600 x 1000	600 x1 300	600 x 1500
Stainless steel wire	ETF70	ETF100	ETF130	ETF150
Stainless steel	ETT70	ETT100	ETT130	ETT150
SHELVES	450 x 700	450 x 1000	450 x1 300	450 x 1500
Stainless steel wire	ETF70.450	ETF100.450	ETF130.450	ETF150.450
Stainless steel	ETT70.450	ETT100.450	ETT130.450	ETT150.450

Liaison d'angle

They allow the mounting of shelves in herring-boat pattern. The link between the 2 shelves requires a pair of connections.

For example, on an angle with 5 shelf levels, you need 5 pairs of angle links.

Ref. 100967 Angle link – pair



FIXED SHELVES WITH EXTENSION

On basic fixed shelving you can add one or more extensions based on your needs, size of your rooms, walking and safety areas.

Basic shelving :

Basic shelves sizes	600x700	600x1000	600x1300	600x1500
Height 1700mm	RF17070	RF17010	RF17013	RF17015
Height 2000mm	RF20070	RF20010	RF20013	RF20015
Basic shelves sizes	450x700	450x1000	450x1300	450x1500
Height 1700mm	RF17070.450	RF17010.450	RF17013.450	RF17015.450
Height 2000mm	RF20070.450	RF20010.450	RF20013.450	RF20015.450

Consisting of : 2 ladders
4 wire shelves (height 1700mm), 5 wire shelves (height 2000mm)
2 support struts

Alignment extension :

Linear extension module	600x700	600x1000	600x1300	600x1500
Height 1700mm	AL17070	AL17010	AL17013	AL17015
Height 2000mm	AL20070	AL20010	AL20013	AL20015
Linear extension module	450x700	450x1000	450x1300	450x1500
Height 1700mm	AL17070.450	AL17010.450	AL17013.450	AL17015.450
Height 2000mm	AL20070.450	AL20010.450	AL20013.450	AL20015.450

Consisting of : 1 ladder
4 wire shelves (height 1700mm), 5 wire shelves (height 2000mm)
2 support struts

Right angle extension :

Herring-boat pattern extension module	600x700	600x1000	600x1300	600x1500
Height 1700mm	EP17070	EP17010	EP17013	EP17015
Height 2000mm	EP20070	EP20010	EP20013	EP20015
Herring-boat pattern extension module	450x700	450x1000	450x1300	450x1500
Height 1700mm	EP17070.450	EP17010.450	EP17013.450	EP17015.450
Height 2000mm	EP20070.450	EP20010.450	EP20013.450	EP20015.450

Consisting of : 1 ladder
4 wire shelves (height 1700mm), 5 wire shelves (height 2000mm)
Pairs of connections

MOBILE SHELVING

Mobile shelving size	600x700	600x1000	600x1300	600x1500
Overall 1110mm	RR11070	RR11010	RR11013	RR11015
Overall 1855mm	RR18570	RR18510	RR18513	RR18515
Mobile shelving size	450x700	450x1000	450x1300	450x1500
Overall 1110mm	RR11070	RR11010	RR11013	RR11015
Overall 1855mm	RR18570	RR18510	RR18513	RR18515

Consisting of : 2 ladders on rollers
 4 support struts

NB : Shelves are sold separately

Shelves sold separately

Designed to be positioned directly on mobile shelving.

SHELVES	600x700	600x1000	600x1300	600x1500
18/10 Stainless steel wire	ETF70	ETF100	ETF130	ETF150
18/10 Stainless Steel	ETT70	ETT100	ETT130	ETT150
SHELVES	450x700	450x1000	450x1300	450x1500
18/10 Stainless steel wire	ETF70.450	ETF100.450	ETF130.450	ETF150.450
18/10 Stainless Steel	ETT70.450	ETT100.450	ETT130.450	ETT150.450



OPERATING AREA TROLLEYS AND ACCESSORIES





Agencinox 60

Technical characteristics

For convenient, modular and mobile storage of your consumables.

In the operating area, these trolleys allow your operating theatres to be supplied daily. Refer to the chapter « STERILIZATION ACCESSORIES » to customize the layout of these trolleys based on your surgical specialty and your organization.

Made of square pipe with a 25X25mm section, 18/10 stainless steel, fully welded.

4 antistatic castor wheels with a diameter of 125mm including 2 with brakes.

SINGLE-SIDED TROLLEY

Reference	Outside dimensions	Inter-hook spacing	Accessory levels
ST83	L.555 x P.450 x H.1396mm	320mm	3 accessoires
ST84	L.555 x P.450 x H.1716mm	320mm	4 accessoires



DOUBLE-SIDED TROLLEYS

Reference	Outside dimensions	Inter-Hook spacing	Accessory levels
ST94	L.555 x P.560 x H.1076mm	320mm	4 accessoires
ST96	L.555 x P.560 x H.1396mm	320mm	6 accessoires
ST98	L.555 x P.560 x H.1716mm	320mm	8 accessoires



DOUBLE-SIDED TROLLEYS WITH CENTRAL COMPARTMENT

Trolley equipped with 4 levels of accessories on each side and a central space designed for baskets and containers with maximum dimensions of 600x300mm.
16 hooks allow an adjustment at every 148mm of the central shelves.
4 castor wheels with a diameter of 125mm including 2 with brakes.

Central shelves ref. ST114 sold separately

Reference	Outside dimensions	Inter-hook spacing	Accessory levels
ST117	L.555 x P.905 x H.1760mm	320mm	8 accessories
ST114	540 x 369mm	Central shelf	



NOTE : All accessories presented in chapter « STERILIZATION ACCESSORIES » are compatible with this trolley

Technical characteristics

In 18/10 stainless steel pipe, with a 25x25mm section and end fittings.

The wall brackets are equipped with welded stainless steel hooks for baskets and accessories with a height of 130 mm to 400 mm based on models.

The brackets are mounted on the wall using stainless steel struts, leaving a space between the wall and the support that facilitates cleaning.

WALL BRACKETS			
Reference	Outside Dimensions	Inter-hook spacing	Accessory level
RM9.2	900mm	589mm	2
RM9.4	900mm	246mm	4
RM13.3	1 300mm	494mm	3
RM13.6	1 300mm	228mm	6
RM16.4	1 600mm	422mm	4
RM16.7	1 600mm	235mm	7
RM17.5	1 700mm	359mm	5

Sold in pairs



Note : All accessories in the chapter « STERILIZATION ACCESSORIES » are compatible with these brackets.

Stainless steel shelves



Multipurpose stainless steel wire shelf, mesh 30x30 with front rim.

ST103 600x300x H150 mm

ST109 600x180x H150 mm



Multipurpose stainless steel wire shelf, stainless steel sheet bottom and front rim.

ST120 dim. 600 x 300 x Ht 150 mm



Multipurpose stainless steel wire shelf, stainless steel sheet bottom, front rim and 30x30mm back mesh.

ST123 dim. 600 x 300 x Ht 260 mm



Double stainless steel wire multipurpose shelf, stainless steel sheet bottom, front and sides rim.

Inside dimensions 590x165mm.

ST124 dim. 600 x 200 x Ht 300 mm



Stainless steel console table with front rim

ST134 600 x 400 mm



Writing desk with 18/10 stainless steel rack

ST111 600 x 350 x Ht 200 mm




18/10 stainless steel wire probe rack,
1 compartment .

507

Dimensions

Width 160 x Height 395 x Depth 70 mm



18/10 stainless steel wire probe rack,
3 compartments 120mm.

507T

Dimensions :

Width 380 x 510 x Depth 70 mm



18/10 stainless steel wire glove rack with 6 compartments

ST62



8 hook holder for probes or anesthesia pipes made of 18/10 stainless steel to be mounted in front or on the side of baskets.

ST132



18/10 stainless steel shoes holder for 6 suspended pairs.

ST61



18/10 stainless steel shoes shelf for 6 pairs.

ST60



18/10 stainless steel wire
decontamination tank or container
holder 530x325 with lid rack storage

ST122



18/10 stainless steel wire container
holder.


ST125



18/10 stainless steel wire container
holder


ST110 : 10 boxes on 2 levels

ST112 : 11 boxes on 4 levels



18/10 stainless steel tubing dispenser.
Dimensions 600x180x150mm for rolls
of 180mm of diameter

ST104



Block with 5 swinging boxes (st.st holder), inside dimensions of a box : 100x100x155mm

ST118

Block with 6 swinging boxes (st.st holder), inside dimensions of a box : 70x100x155mm

ST119



Holder for spout pans with a capacity of 1l, 3l, 7l, 10l. On a single level, you may have :

> 5 1-liter pans

> 4 3-liter pans

> 3 7-liter pans

ST139



1-liter spout pan. Overall dimensions : width 102, depth 157, height 73mm
ST141 Blue.

3 -liter spout pan. Overall dimensions : width 145, depth 280, height 180mm
ST142 Blue.

7 -liter spout pan. Overall dimensions : width 210, depth 340, height 145mm
ST143 Blue.

10 -liter spout pan. Overall dimensions : width 204, depth 340, height 198mm
ST144 Blue.

BASKETS

Reference	Name and outside dimensions
Non-notched baskets	
ST101	Stainless steel wire basket 600x300x300mm
ST102	Stainless steel wire basket 600x300x150mm
ST105	Stainless steel wire basket 600x400x300mm
Notched baskets	
ST106	Stainless steel wire basket 600x300x300mm
ST107	Stainless steel wire basket 600x300x150mm
ST108	Stainless steel wire basket 600x400x300mm

18/10 Stainless steel wire, sturdy.
Rounded finish to prevent tearing of pockets or papers
Side handgrip for easy handling.
Stackable baskets
Standardized models



ST51

BASKET ACCESSORIES

References	Outside dimensions
WIDTH DIVIDERS	
ST126	Divider 300 x 150 mm
ST127	Divider 300 x 300 mm
ST129	Divider 400 x 300 mm
LENTH DIVIDERS	
ST131	Divider 600 x 300mm
ST135	Divider 600 x 150 mm



Width divider

DECONTAMINATION TROLLEYS AND ACCESSORIES





ST153

PRION CART

ST153

Totally made from stainless steel 18/10 aisi 304 grade. Frame is made from square stainless steel tube 25x25mm, and argon welded. Fitted with 2 stainless steel basins without any retention, with lid and draincock lower shelf is removable covering the frame and fixing parts.

4 castors Ø100mm including 2 with brakes.

PREDESINFECTION TROLLEY ST154

Totally made from stainless steel 18/10 aisi 304 grade. Frame is made from square stainless steel tube 25x25mm, and argon welded.

Fitted with 40L high edges container with lid and draincock.

4 castors Ø100mm including 2 with brakes.

Overall dimensions 845 x 475 x Ht1020mm



ST154

CONTAINER TROLLEYS

Totally made from stainless steel 18/10 aisi 304 grade. Round tube frame Ø25mm and argon welded. Removable stainless steel lower shelf. Fitted with 60L polypropylene container with lid.

4 castors Ø100mm including 2 with brakes.

Option : Draincock to be specify in the order

Designation	Ref.
Container trolley Dimensions 750x450xHt940	ST152
Container trolley top moved to the front for easy draining. Dimensions 720x450xHt940	ST155
Draincock	ST838



ST152



ST155

This trolley line was especially designed to facilitate transportation of decontamination pans and instrument containers from operating area to sterilization area..

In order to meet your needs, based on your organization and restrictions (lifting platform, storage space, etc...) we provide several types of equipment.

PLATFORM RACK TROLLEY

Designed to accommodate 530x325mm decontamination pans, this trolley made of 18/10 stainless steel is equipped with an interesting modularity. The 3 pairs of rails, welded for better stability of the unit, allow you to position :

- 3 pans, height 100mm
- 3 pans, height 150mm
- 2 pans, height 200mm

In addition, you may position stainless steel trays on the rails to transport boxes or instrumentation containers or even decontamination pans with dimensions that are less than 530x325mm.

In order to provide optimal mobility, this trolley is equipped with 4 wheels with a diameter of 125mm including 2 with brakes. To clean it in a washer, we suggest the option with special wheels for washer.

Reference	Name	Outside dimensions (mm)
ST81	Decontamination trolley	570 x 590x Ht 1000-- mm
ST80	Removable stainless steel tray	



POLYCARBONATE

Technical characteristics

Polycarbonate is a thick translucent smoke-coloured recycling plastic material, withstanding temperatures from -40° to +100°C. As a result, this material is compatible with cleaning in a washer.

Highly resistant to shocks and rigid due to exceptional thickness, this material is in compliance with all food and hospital hygiene standards, is a NF approved material and meets the criteria related to European standard EN631.1

Polycarbonate Decontamination pans

The pans below have excellent ergonomics : easy to stack and light, these pans are all the more easy to use as they are graduated and designed for quick drying and optimal draining.



Reference	Outside Dimensions	Height (mm)	Capacity
ST822	530 X 325	100	13 L
ST821		150	19,5 L
ST823		200	25,7 L
ST820		Lid with handle	
ST824		Sealed lid	

POLYPROPYLENE

Technical characteristics

Polypropylene is a thick opaque plastic material, withstanding temperatures from -30°C to +100°C. Although this material is compatible with cleaning in a washer, it cannot be passed through an autoclave at temperatures over 100°C.

More flexible than polycarbonate, polypropylene is more resistant to shock.



Polypropylene decontamination pans

Reference	Outside dimensions	Height (mm)	Capacity
ST834	530 x 325	100	13 L
ST835		150	20,5 L
ST836		200	27,2 L
ST837		Polycarbonate transparent lid	
ST838		Drainage Kit	
ST833		Disinfection Kit	

Drainage Kit

Because we want to make decontamination easier for you, we provide on these pans as an option a drainage kit : tap, nut and connection.



Disinfection Kit

A polypropylene pan of 530x325mm and height 150mm have a capacity of 20,5 liters.
A PVC thermoformed lid.
A draining rack.



Polypropylene pans for endoscope

Opaque pans for endoscope dipping. Sealed lid. Handle facilitating handgrip.



Reference	Outside dimensions	Height (mm)	Capacity
ST839	560 x 355	150	32.2 L
ST840		Lid	

DURADEX : Duradex decontamination pans

Technical characteristics

Duradex is a resin designed to withstand stains and scratches. Its durability is guaranteed by increased resistance to temperature variations as well as chemical treatments.

Duradex pans are compatible with washer and even autoclave cleaning, and because they are resistant to temperatures between -40°C and +200°C, they can be subjected to a Prion cycle.

Transparent, amber colored, these pans are perfectly identifiable in order to facilitate compliance with your supply items.

The Duradex lids are subjected to the same treatments as their pans.



Reference	Outside dimensions	Height (mm)	Capacity
ST826	530 x 325	100	13 L
ST827		150	19,5 L
ST828		Duradex Transparent Lid	

STAINLESS STEEL : Stainless steel Decontamination

Technical characteristics :

In order to guarantee the long-life of our line of 18/10 stainless steel pans, the corners and edges of our pans are rigidified by stamping. As a result their resistance to impacts and falls is increased and the risks of distorsion are reduced.

Reference	Outside Dimensions	Height (mm)	Capacity
ST841	530 x 325	100	13,5 L
ST842		150	20 L
ST843		200	28 L
ST844		Lid with handle	
ST845		Sealed lid	

All stainless steel references

PACKING TABLES





Agencinox 78


PACKING TABLES

Sturdy structure in 18/10 stainless steel, fully welded.

> Anti-reflective scratch-brushed stainless steel tray with a thickness of 1.5 mm with 35mm rim, at an angle of 45° for optimal cleaning and hygiene.

> Compact tray class M1 lab quality Tacheeclean (4 colors available : Aquamarine, Oslo blue, Nectarine, Grapefruit)

4 castor wheels of 125mm diameter, 2 with brakes.



TABLES	1200x700mm	1500x700mm	1800x700mm
Stainless Steel tray	STC120	STC150	STC180
Compact tray	STC120C	STC150C	STC180C



OPTIONS

Cross-members

To hang packing accessories.

Made of bent square pipe 25x25mm, 18/10 stainless steel

Equipped single-sided (2 levels of accessories) or double-sided (4 levels of accessories)

Two cross-members may be adapted side by side

Cross-member height : 1,487 mm, length 555 mm.

Cross-members	Reference
Single-sided 2 levels	STP21
Double-sided 4 levels	STP22

Pan support blocks

Mounted on the working surface on the right or on the left side, our pan holder is adjustable and can be equipped with 2 polycarbonate pans 600x400x200 mm or 3 pans 600x400x100 mm

Reference	Name
ST148	Pan support block
6200	Pan 600 x 400 x 200 mm
6202	Pan length divider 6200
6203	Pan width divider 6200
6100	Pan 600x 400 x 100 mm
6102	Pan length divider 6100
6103	Pan width divider 6100

PAPER SUPPORT TROLLEY

ST99 Paper support trolley Made of stainless steel

4 bars accommodating all paper formats. Paper support system.


Lower stainless steel wire grid. 4 antistatic wheels dia. 100mm including 2 with brakes.

Overall dimensions : 1,330x650mm. Height 970mm



TRANSFER AND EMERGENCY STRETCHERS WHEEL CHAIR

AXIOME
Aviso
ABSOLUTE
AMBULYS



Comfort

Safety

CE

Reliability

Variable height



- Height adjustable by hydraulic spring, left and right operated
- Removable mattress. Flame retardant M1 (83% PVC / 17% PES). Resin hard top
- Collapsible side rail
- Adjustable backrest (from 0 to 90°). Double gaz spring assisted.
- Telescopic and collapsible serum rod – 2 hooks.
- Accessories holders and bumpers.
- 4 wheels 200 mm diam. Optional 5 th wheel for better manoeuvrability
- 2 retractables handles
- Epoxy (grey 7035 and 9007). Corrosion resistant – tropicalized
- Optional trendelenburg position by double gaz spring .
“Easy lift system”

AXIOME

Variable height	Trendelenburg
Width 500mm	Width 500mm
CB501	CB502
CB501MC	CB502MC
CB501 5R	CB502 5R
CB501MC5R	CB502MC5R
Width 600mm	Largeur 600mm
CB601	CB602
CB601MC	CB602MC
CB601 5R	CB602 5R
CB601MC5R	CB602MC5R
Width 700mm	Width 700mm
CB701	CB702
CB701MC	CB702MC
CB701 5R	CB702 5R
CB701MC5R	CB702MC5R
Width 800mm	Width 800mm
CB801	CB802
CB801MC	CB802MC
CB801 5R	CB802 5R
CB801MC5R	CB802MC5R

Trendelenburg




Knee section

Knee section
Mattress thickness 120mm
Width 600mm
CB602P
CB602P 5R
Width 700mm
CB702P
CB702P 5R


Accessories

- **Ref. CBSO** Oxygen holder
- **Ref. CBPA** Basket
- **Ref. CBSS** Monitoring holder with straps
- **Ref. CBSR** Paper roll holder
- **Ref. CBSU** Urinal bag holder
- **Ref. CBTSA** Additionnal telescopic serum rod
- **Ref. CBTSR2** Additionnal telescopic foldable serum rod (*head end only*)
- **Ref. CBTR** Food tray
- **Ref. CBPR** File holder
- **Ref. CBK7** Radio holder including the full length Translatory motion
- **Ref. CBSM** Over bed radiotransluscent mattress
- **Ref. CBARE** Shoulders restraining strap
- **Ref. CBARM** Ankle restraining strap
- **Ref. CBARP** Wrist restraining strap
- **Ref. CBPB** Side rail protection
- **Ref. CBT** Ophtalmic surgery headrest
- **Ref. CBPP** Additional pushing handle



Technical specifications

Mattress width	Mattress length	Mattress thickness	Mattress Thickness (knee sect.)	Overall width	Overall length	Low position	Hight position	Weight without accessories	Load capacity
500mm	1860mm	80mm		620mm	2040mm	580mm	890mm	70kg	250kg
600mm	1860mm	80mm	120mm	720mm	2040mm	580mm	890mm	80kg	250kg
700mm	1860mm	100mm	120mm	820mm	2040mm	580mm	890mm	100kg	250kg
800mm	1860mm	100mm		920mm	2040mm	580mm	890mm	110kg	300kg



Model presented : Aviso Ref. CBL601MC with CBLPPE



Model presented : Aviso Ref. CBL6025R with CBLTS and CBLPPE with trendelenburg

Configuration and specifications

- Height adjustable by hydraulic spring.
- Trendelenburg version available.
- Removable mattress, flame retardant M1, 35 kg/m³ housse 83% PVC / 17% PES, resin hard top.
- Collapsible side rail.
- Adjustable backrest (from 0 to 90°). Double gaz spring assisted (easy lift).
- Accessories holders and bumpers.
- 4 wheels Ø150 mm. Central braking system left/right side operated. (optional wheels Ø200mm)
- Option : 5th wheel for better manoeuvrability
- Epoxy coated frame (grey 7035 and 9007). Corrosion resistant - tropicalized.
- 200 kg capacity.

Transfer mattress available, PVC coated coton cover.

Optional accessories

- | | |
|--|---------------------|
| • Basket / Shelf | Ref. CBLPA |
| • 2 retractables handles
(head or foot end) | Ref. CBLPPE |
| • 2 fixed handles
(head or foot end) | Ref. CBLPPF |
| • Removable serum rod | Ref. CBLTS |
| • Telescopic and foldable serum rod
(head end only) | Ref. CBLTSR2 |
| • Removable telescopic serum rod | Ref. CBLTSA |
| • Holder for one O2 | Ref. CBLSO |
| • Urinal bag holder | Ref. CBLSU |
| • Chart/radio holder | Ref. CBLPR |
| • Paper roll holder | Ref. CBLSR |
| • 2 aluminium rail | Ref. CBLRA |

Dimensions

Height	Trendelenburg
Width 600mm	Width 600mm
CBL601	CBL602
CBL601MC	CBL602MC
CBL601 5R	CBL602 5R
CBL601MC5R	CBL602MC5R
Width 700mm	Width 700mm
CBL701	CBL702
CBL701MC	CBL702MC
CBL701 5R	CBL702 5R
CBL701MC5R	CBL702MC5R

Technical specifications

Mattress width	Mattress length	Mattress thickness	Overall width	Overall length	Low position	High position	Weight without accessories
600mm	1860mm	80mm	717mm	2040mm	525mm	917mm	80kg
700mm	1860mm	100mm	817mm	2040mm	525mm	917mm	95kg



Comfortable pushing handle

Electrical backrest


Control pad 4 adjustments



Shelf utility tray

Ø 100 wheel with brake

Electrical leg-rest



FA.GF
Grey / framboise



FA.LP
Lavande / prune



FA.GO
Grey / orange



FA.AM
Anis / brown

Configurations and specifications

Patient chair, electrical with 4 « LINAK » actuators allowing height adjustment from 560 mm to 920 mm, Trendelenburg / reverse 20°/18°, backrest and legrest from 90° to 180°. All adjustments are done by hand control pad.

Features:

- Basket for accessories.
- 4 x 100mm wheels including 2 with brake. (non marking).
- Epoxy grey coated 7035 RAL , frame finishing, anti rust treatment, corrosion resist tropicalized.
- Max. load 200 kg.
- Saddle : choice of 4 bi-colors upholstery.
- Plastic cover for leg-rest

3 parts couch with additional memory foam, 100mm thick 35 kg/m³, STAMSKIN® coating upholstery: Abrasive resistant testing Martindale 120 000 cycles.

Clean regularly with a damp soapy cloth. Do not use strong detergents, solvents cleaners, polish, wax or patent chemicals: harmful for the material. Some stains like ball pen, food coloring should be cleaned off immediately to avoid any permanent absorption into the upholstery.

■ Technical features

STAMSKIN
TOP



Confort and flexibility



Tear resistance



Long life



Waterproof



Easy care



Fire proof norm
M2 (NMP 92.507)



Fungicide
treatment



Recyclable
textile 100%

Optionnal accessories (assembly plant)

- | | |
|---|------------------------|
| • Revolving arm-rest | Ref. FAAR |
| • Adjustable blood taking arm-rest on ball | Ref. FAABR |
| • 2 stainless steel rail 25x10mm | Ref. FARA |
| • Paper roll holder | Ref. FASR |
| • Urinal bag holder | Ref. FASU |
| • Battery | Ref. FABAT |
| • Stainless steel foot-rest | Ref. FARP |
| • 4X125mm directional and central braking
Wheels | Ref. +color ref |
| • Additional Brake on front wheels | Ref. FAPF |

Dimensions

Couch :

- Length 1900mm
- Width 550mm

Ext. dimensions:

- Length 1900mm
- Width 820mm
(arm-rest version)

Height:

- Position low 560mm
- Position high 920mm

Other optionnal accessories

Gross weight : 80 kg

- | | |
|-------------------------------|------------------|
| • Armrest clamp Ø 16mm | Ref. FAEF |
| • Multi use shelf | Ref. FAPC |
| • Serum rod 2 hooks | Ref. FATS |
| • Flexible control pad holder | Ref. FAST |
| • Adjustable head pillow | Ref. FART |
| • LED reading light | Ref. FALL |

AMBULYS

MULTI PURPOSE CHAIR
Mattress size 600 or 700mm



Pushing handle,
 Seat or flat
 position



Safety side rail and
 extensible 1000 to
 1150mm (option)

Control pad 4
 adjustments

Electrical
 back-rest



Emergency
 position control

Electrical
 leg-rest

Ø125 wheel or Ø150
 (with brake)



FHJ.GF
 Grey / framboise



FHJ.LP
 Lavande / prune



FHJ.GO
 Grey / orange



FHJ.AM
 Anis / brown

Configuration and specifications

- 4 electrical adjustments with control pad : Height, back-rest, leg-rest, trendelenburg, reverse trendelenburg
- Pushing handle seat and/or flat position.
- Removable safety side rail length 1000mm or adjustable side rail from (1000 to 1150mm), in accordance with the standard bed-ins, useful for adults and children.
- Emergency position control.
- Accessories holder.
- 4 wheels with brake Ø 125mm 2 side control (no marking).
- Epoxy grey coated 7035, corrosion resist tropicalized.
- 200 kg loading capacity.
- Saddle : choice of 4 bi-colors upholstery.

Saddlery with memory foam depth 100mm 35 kg/m³, coating STAMSKIN® : Abrasive resist testing Martindale 120 000 cycles. Plastic cover for leg-rest

Clean regularly with a damp soapy cloth. Do not use strong detergents, solvents cleaners, polish, wax or patent chemicals: harmful for the material. Some stains like ball pen, food colorings should be cleaned off immediately to avoid any permanent absorption into the vinyl.

Leg-rest protection.

■ Technical features

STAMSKIN
TOP



Confort and flexibility



Tear resistance



Long life



Waterproof



Easy care



Fire proof norm
M2 (NMP 92.507)



Fungicide
treatment



Recyclable
textile 100%

Optionnal accessories (assembly plant)

- 2 rails inox 25x10mm **Ref. FHJRA**
- 4 wheels Ø150mm central braking **Ref. FHJR150**
- Urinal bag holder **Ref. FHJSU**
- Battery **Ref. FHJBAT**
- Stainless steel foot-rest **Ref. FHJCPF**
- Stainless steel adjustable foot-rest **Ref. FHJCPR**

Specifications / Dimensions


	FHJ600	FHJ700
Mattress		
■ Length	2100	2100
■ Width	600	700
Overall		
■ Length	2190	2190
■ Width	780	880
height		
■ Low position	540	540
■ High position	905	905
Weight	90kg	96kg

Other optionnal accessories

- Clamp Ø 16mm **Ref. FHJEF**
- Multi use shelf **Ref. FHJPC**
- Serum rod 2 hooks **Ref. FHJTS**
- Flexible control pad holder **Ref. FHJST**
- Adjustable head pillow **Ref. FHJRT**
- LED reading light **Ref. FHJLL**
- 2 armrest mounted on ball **Ref. FHJAB**
- Brake on front wheel side **Ref. FHJPF**

ABSOLUTE **AMBULYS**

Comfortable - Safety - Ergonomic - Handling hability - Reliability - Sturdiness



WHEEL CHAIR OPTIONS





CARTS and CABINETS MEDICATION and STORAGE



Ref. DMCH011

Cart / Trolley supplies



Medication cart :

The range of medication trolleys provides the best maneuverability and optimum impact resistance ; prevent from chemicals attacks thanks to the selection of raw material used.

The single or double version have been designed to optimize the space and time (for medication tasks) due to the small boundary dimensions and high capacity with a range from 16 to 36 weekly trays and number of baskets 400x300mm.

Dimensions :

Length: 720mm
Depth : 520mm
Height : 1270mm

Specifications :

- Top shelf rack from ABS
- Resin inside wall
- Rails made from ABS
- Castors Ø 125mm with 2 brakes
- Closing with PVC rolling shutter
- Safety key locking system
- Anodised aluminium frame
- Adjustable handle and accessories
- Optional electronic locking key (ref. DMAC006)



Ref. DMCH012

Trolley for 32 pill dispenser trays

Dimensions :

Length : 1145mm
Depth : 530mm
Height : 1270mm

Ref. DMCH013

Trolley for 36 pill dispenser trays

Dimensions :

Length : 1145mm
Depth : 530mm
Height : 1400mm

Presentation with accessories.

Ref. DMAR003

Storage cabinet :

The cabinet allows to fit easily the hospital pharmacies in accordance with hygiene rules.

Easy to manage the daily procurement for medication in the different hospital departments thanks to the compatibility with the carts range (ref. DMCH011, DMCH012 et DMCH013 customized for pillbox trays ref. DMPI015) or baskets with dividers 400x300mm(DMPI058 and DMPI059)

Dimensions :

Length : 540mm
Depth 380mm
Height : 1950mm

Specifications :

- Epoxy coated steel frame (blue or white)
- PVC rolling shutter
- Key locking safety system
- ABS inside wall easy cleaning
- Capacity of 30 pillbox trays

Optional electronic locking key (ref. DMAC006)

ACCESSORIES

Daily pill dispenser :

Essential and convenient for the daily medication.

Set of 7 pill boxes fitting in a tray for a week including 4 adjustable and removable day dividers, identification stickers.

The number of accessories allow to customize the medication to the hospital practice.

Cleaning temperature up to 50°C



DMPI015

Dimensions :
Length : 400mm
Depth : 300mm
Height : 47mm



DMPI058

ABS basket
400x300x100mm



DMPI020

Tray for pill dispenser



DMPI018

Day pill dispenser



DMPI016

Day stickers



DMPI019

Sticker holder





Ref. DMGH002

Dimensions :
Length : 880mm
Depth : 530mm
Height : 1120mm
Unladen weight : 34kg



Ref. DMGH014

Dimensions :
Length : 880mm
Depth : 530mm
Height : 1270mm
Unladen weight : 36kg

Carts and mobile storage cabinet ISO 600x400 :

The modular range of carts and cabinets gives to the users many possibilities for medication and consumables storage in the hospitals,
The multiple configurations are made by the wide range of accessories, like baskets or dispenser trays allowing to comply with the different medication and distribution process.
Sturdy construction made in accordance with ergonomic constraints and users requirements.



Ref. DMAR001

Specifications :

- Top shelf made from ABS
 - Resin inside wall
 - Rails made from ABS
 - Castors Ø125mm, 2 with brakes
 - Closing with PVC rolling shutter
 - Locking by key with safety system
 - Anodised aluminium frame
 - Adjustable handle and accessories
- Optional electronic locking key (ref. DMAC006)

Dimensions :
Length : 880mm
Depth : 505mm
Height : 1850mm
Unladen weight : 58kg

Storage cabinet ISO 600 X 400 :

The cabinet is made for pharmacies storage area and fits with all ISO 600 x 400 (pill dispenser trays, baskets, ISO trays...) safety drug cabinet included lockable by key.

- Blue epoxy coated frame
- PVC rolling door
- Key locking safety system for main door and drugs
- Inside wall and rails made from ABS
- 22 pill dispenser trays capacity (DMPI021)
- Electronic locking key in option (DMAC006)


Dimensions :
Length : 750mm
Depth 500mm
Height: 1960mm



Ref. DMAR004

Presentation with accessories.

ACCESSORIES ISO 600 x 400



DMPI057
Tray



DMPI021
10 pillbox tray



DMPI024
Pillbox
400x92x82



DMPI025
Pillbox
400x132x100




Ref. DMPI050
ABS basket
Length : 600mm
Depth 400mm
Height : 50mm



Ref. DMPI052
ABS basket
Length : 600mm
Depth 400mm
Height : 100mm



Ref. DMPI054
ABS basket
Length : 600mm
Depth 400mm
Height : 200mm



DMAC018
Accessories holder



DMAC002
Adjustable handle



DMAC004
Glove box holder



DMAC007
Telescopic shelf



DMAC003
Waste container



DMAC005 Sharps
container holder



DMAC009
Basket for bottles

DMAC012
Serum rod
stainless steel



DMAC006
Electronic locking key



DMAC010
Laptop holder



DMAC016
Retractable side shelf



DMAC011
Overbridge unit for tilt out bin

Ref DMCH020

Anesthesia cart :

Ergonomic and hygienic with adjustable handle and accessories. Easy for cleaning and disinfection. Top shelf made from ABS. 6 drawers opening at 110%, fitting ISO baskets easily removable mounted on telescopic rails. Ø125 wheels, with 2 brakes. Central locking by key.

Accessories :

- Overbridge unit with storage boxes
- Waste bin easy to remove and clean
- Needle collector holder
- Glove box holder
- Drawers dividers

Unladen weight : 36kg

Dimensions :
Length :880mm
Depth 525mm
Height: 1072mm



Ref DMCH021

Emergency cart :

Standard configuration based on the anesthesia cart except the seal locking drawers.

Accessories :

- Heart massage board
- IV pole
- Defibrillator holder
- Oxygen bottle holder
- Waste bin with holder
- Needle collector holder
- Glove box holder
- Drawer dividers
- Seal locking system

Dimensions :
Length :880mm
Depth 525mm
Height: 1072mm



Presentation with accessories.

Ref. DMGH001

Multifunction trolley

Specifications :

- ABS frame (blue and white)
- 9 drawers with removable dividers
- Key locking safety system
- Castors Ø 120mm with 2 front brakes
- Easy cleaning

Dimensions :

Length : 850mm
Depth 640mm
Height: 1150mm

Waste container and accessories are available on page 5.

Presentation with accessories.



Ref. 980

Storage trolley

Specifications :

- Stainless steel frame
- Castors Ø 100mm, 2 with brakes, with bumpers
- 18 adjustable positions

Delivered without basket and rails

Réf. 6606 Stainless steel horizontal rails x2 length 600mm

Réf. 6608 Stainless steel angle rails x2

Dimensions :

Length: 455mm
Depth 640mm
Height: 1670mm



Ref. 885

Transporting trolley :

Specifications :

- Stainless steel frame
- Castors Ø100mm with 2 brakes, with bumpers

Delivered without basket and rails

Réf. 6606 Stainless steel horizontal rails x2 length 600mm

Réf. 885P Stainless steel top tray

Dimensions :

Length: 885mm
Depth 640mm
Height: 1020mm




Model présenté with stainless steel top shelf, stainless steel rails and baskets.

Baskets are available on page 5

AGENCINOX

Your distributor :



© Agencinox 2013

AGENCINOX

Z.I. 89110 Aillant sur Tholon

S.A.S. RC Sens B 353 427 495 – Code APE 3250 A
Export : Tel. 33 3 86 92 39 27 – Fax 33 3 86 91 54 42

contact@agencinox.com

www.agencinox.com

AGENCINOX certified



POLYCERT

ISO 9001 : 2008
ISO 13485 : 2012





**WASH-ARMS AND
WASH-HAND BASINS**

THE XB RANGE



XB 14

Knee control, flow of mixed water controlled by knee lever

Basin :

Width : 405 mm
Projection : 405 mm
(with backboard : 421 mm)
Depth : 120 mm

With backboard equipped with soap dispenser !



XB 15

Knee control, flow of mixed water controlled by knee lever.



XB 24

Leg control, flow of mixed water controlled by leg flap.

With backboard equipped with soap dispenser !



XB 25

Leg control, flow of mixed water controlled by leg flap.



XB 34

Electronic control, flow of mixed water controlled by infrared cell

XB34P

With two Lithium battery 6 Volts. Easy installation

With backboard equipped with soap dispenser !



XB 35

Electronic control, flow of mixed water controlled by infrared cell

XB35P

With two Lithium battery 6 Volts. Easy installation

WASH-ARMS AND WASH-HANDS BASINS

WALL-MOUNTED BACKPIECE

Model on picture : MAX 2S

The whole of this range is designed to provide you with complete equipment in cases of little space. Its backpiece facilitating the installation of accessories is fitted with a tight seal in view of simple and rapid installation.



MAX 2S

Wash-hand basin 18/10 stainless steel. Knee control of water, towel dispenser and soap dispenser in smoked ABS.

WALL-MOUNTED BACKPIECE WASH-ARM BASIN

Model in picture : ELECTRO MS

The whole of this range is designed for intensive use. Through the space between the height of the water swanneck and the large capacity tank, it is possible to very easily wash the forearms.



ELECTRO MS

Wash-arm basin 18/10 stainless steel. Electronic tap, towel dispenser and soap dispenser in smoked ABS.

LMG 2S

Wash-arm basin 18/10 stainless steel. Knee control of water, towel dispenser and soap dispenser in smoked ABS.

THE BASIC BASINS



MAX WALL-MOUNTED BACKPIECE WASH-HAND BASIN

18/10 Stainless steel, medium capacity tank, rounded front corners inclined counterplies for easy cleaning.

Width : 405 mm
Projection : 425 mm
Backpiece height : 550 mm
Tank depth : 120 mm



XB WALL-MOUNTED WASH-HAND BASIN

18/10 Stainless steel, medium capacity tank.

Width : 405 mm
Projection : 405 mm
Tank depth : 120 mm

LBM WALL-MOUNTED BACKPIECE WASH-HAND BASIN

18/10 Stainless steel, high capacity tank, rounded front corners inclined counterplies for easy cleaning..

Width : 500 mm
Projection : 565 mm
Backpiece height : 450 mm
Tank depth : 175 mm



LIZA WASH-HAND ROUNDED BASIN

18/10 Stainless steel, medium capacity tank, inclined counterplies for easy cleaning.



THE ROUND WASH-HAND BASIN MADE IN 18/10 STAINLESS STEEL IS DESIGNED TO BE EASILY FITTED IN CASE OF SMALL PLACE.



LIZA V

Leg control, flow of mixed water controlled by leg flap.

BASIN :
Width : 405 mm
Projection : 466 mm
Depth : 120 mm

**ELECTRONIC
CONTROL ON
BATTERY
LITHIUM 6 V:
EASY TO SET !**



LIZA E

Electronic control, flow of mixed water controlled by infrared cell. Easy to set.



ACCESSORIES FOR LIZA



THE LIZA RANGE COULD BE EQUIPED WITH TOWEL AND SOAP DISPENSER



Réf. DSM

18/10 Stainless steel back board with soap dispenser 0,5 L smoked ABS.

Réf. DSMP

18/10 Stainless steel back board with soap dispenser 0,5 L and single use hand towels dispenser. Smoked ABS



Doc. N°DC658AG