


ALLERGEN DETECTION KITS

4 STEPS WORKFLOW



BIOPFS-0049 • Real Time Detection Kit Celery

BIOPFS-sp0049 • Real Time Detection Kit Celery SAMPLE

Biopremier food safe meat DNA detection kits are based on 5' nuclease real time PCR reactions to amplify unique genomic sequences and enable the simple and reliable detection of animal species. Also suitable to be used with animal feedstuff and other samples in which animal species detection is needed.

- ✓ **Compatible with instruments working in FAM and ROX channels**
- ✓ **Detects low amounts of target DNA**
- ✓ **Ready to use format**
- ✓ **100% specificity**
- ✓ **100 RXN**

| DESCRIPTION | KIT VOLUME | SAMPLE VOLUME |
|----------------------|------------|---------------|
| Master Mix XIII | 2 x 840 µl | 168 µl |
| 10x Celery Assay Mix | 210 µl | 21 µl |
| Positive Control | 70 µl | 10 µl |
| Negative Control | 70 µl | 10 µl |

For development pipeline contact sales.support@biopremier.com
BPMR is certified ISO 9001:2015