

Stomacher® 80 Biomaster Bags

The Stomacher® Bag Range has been designed to fulfil all applications from blending, straining and storage including blender bags designed specifically for life science applications.

Enhanced **Stomacher® 80 bag** range features closure and strainer bag formats, double blender bags for biosafety and the new microBag for effective homogenisation and sample recovery below 1ml volumes.



Standard Bag
105 x 155mm
5-80ml
200 (sachets of 50)
BA6040



Closure Bag
105 x 155mm
5-80ml
200 (sachets of 10)
BA6040/CLR



Strainer Bag
105 x 155mm
5-80ml
200 (sachets of 10)
BA6040/STR



Strainer Double Bag
105 x 155mm
5-80ml
200 (sachets of 10)
BA6040/STR/DBL



Strainer Closure Bag
105 x 155mm
5-80ml
200 (sachets of 10)
BA6040/CLR/STR



Micro Bag
60 x 185mm
250µl-5ml
50 Bags
BA6040/MICRO