

Tube



Other tube labels are available, please refer to tube sections in this catalogue.



Paper Label

The standard label providing a free text space for writing patient ID & sample information.



Block Label

A formatted paper label with defined areas for writing specific patient ID & sample information.



See Through Label

Product specifications are sprayed directly onto the tube with no provision for hand written information. This system reduces the risk of tubes getting stuck in sample racks when multiple layers of inhouse patient identifier labels are subsequently applied to the tube.



Mylar Label

In the same format as a paper label but with the additional benefit of being see through to allow greater visibility of sample.

Shelf Pack



Case



Product Symbols



The CE mark, signifying compliance with the European IVD MD Directive, 98/79/EC



Catalogue or re-order number



Lot number or batch number



Use by, expires or best before



Use once or do not reuse



Sterilised by moist heat



Sterilised by irradiation



Read instruction leaflet before use



Sterilised by Ethylene Oxide gas



Protect from sunlight (may show temperature range)



Protect from any light source



Fragile



Storage temperature range



This way up



Recycle



Date of manufacture



Serial number



In Vitro Diagnostic Medical Device

K2E

EDTA - Dipotassium salt

K3E

EDTA - tripotassium salt

N2E

EDTA - disodium salt

9NC

Trisodium citrate 9:1

4NC

Trisodium citrate 4:1

FE

Fluoride/EDTA

LH

Lithium Heparin

NH

Sodium Heparin

Z

None (no additive)

SST™ II Advance

Silica (clot activator)/Gel

CTAD

Buffered Citrate 0.109M

RST

Thrombin based medical clotting agent/Gel

ACD-A

ACD Solution A

ACD-B

ACD Solution B

LH PST™ II

Lithium Heparin Gel

NC Ficoll

Sodium Citrate/Ficoll

NH Ficoll

Sodium Heparin/Ficoll